



γ -G-test

cinética AA

Modified Szasz method for the determination of γ -glutamyl transferase in serum or plasma.
Substrate recommended by the IFCC.


SUMMARY

γ -glutamyl transferase (γ -GT) is a membrane enzyme widely distributed in the body. It is primarily located in the kidney, seminal vesicles, pancreas, liver, spleen and brain. Its activity is influenced by any factor altering the cellular membranes of the organs that contain it. In the case of liver disorders, γ -GT generally indicates toxic aggression. However, its determination only has clinical relevance when its values are compared to those of other greater organ-specificity enzymes.

γ -GT determination together with alkaline phosphatase, transaminase and bilirubin, significantly broadens the spectrum for differential diagnosis of primary and secondary liver diseases, being part of the hepatic profile.

PRINCIPLE

γ -glutamyl transferase is a carboxypeptidase that catalyzes the following reaction:



PROVIDED REAGENTS

A. Reagent A: Tris buffer solution containing glycylglycine.

B. Reagent B: L- γ -Glutamyl-3-carboxy-4-nitroanilide solution.

Final concentrations

Tris buffer..... 100 mmol/l; pH 8.5
L- γ -Glutamyl-3-carboxy-4-nitroanilide..... > 2.9 mmol/l
glycylglycine 100 mmol/l

NON-PROVIDED REAGENTS

Saline solution (9 g/l NaCl).

INSTRUCTIONS FOR USE

Provided Reagents: ready to use. They may be used separately or as **Monoreagent**, mixing 4 parts Reagent A with 1 part Reagent B (e.g. 4 ml Reagent A + 1 ml Reagent B).

WARNINGS

Reagents are for "in vitro" diagnostic use.
Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: are stable at 2-10°C until the expiration date stated on the box.

Monoreagent (premixed): stable for 4 weeks at 2-10°C. Protect from sunlight.

INSTABILITY OR DETERIORATION OF REAGENTS

When the spectrophotometer has been set to zero with distilled water, absorbance readings of the Monoreagent higher than 1.300 O.D. indicate reagent deterioration.

SAMPLE

Serum or plasma

a) Collection: obtain sample in the usual way.

b) Additives: when using plasma, collect it with EDTA as anticoagulant.

c) Known interfering substances: no interferences have been observed by bilirubin up to 280 mg/l (28 mg/dl), triglycerides up to 5.4 g/l (540 mg/dl) nor hemoglobin up to 0.39 g/dl (390 mg/dl). See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: γ -GT in serum is stable for 2 weeks at 2-10°C and up to 6 months at -20°C.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer.
- Micropipettes and pipettes for measuring the stated volumes.
- Spectrophotometric square cuvettes.
- Water bath at selected assay temperature.
- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 405 nm
- Reaction temperature: 25, 30 or 37°C
- Reaction time: 3 minutes
- Sample volume: 100 μ l
- Reagent B volume (Substrate): 1 ml
- Final reaction volume: 1.1 ml

PROCEDURE

In a cuvette kept at the selected temperature, place:

Monoreagent 1 ml

Pre-incubate a few minutes. Then add:

Sample 100 μ l

Mix at once and continue incubation immediately and simultaneously start stopwatch. Record absorbance at 1,

2 and 3 minutes. Determine average absorbance change ($\Delta A/\text{min}$) subtracting each reading from the previous one and averaging values. Use this mean for calculations.

CALCULATIONS

γ -Glutamyl Transferase (U/l) = $\Delta A/\text{min} \times 1,158$

SI SYSTEM UNITS CONVERSION

γ -GT (U/l) $\times 0.017 = \gamma$ -GT (ukat/l)

QUALITY CONTROL METHOD

Each time the test is running, analyze two levels of a quality control material (**Standatrol S-E 2 niveles**) with known γ -glutamyl transferase activity.

REFERENCE VALUES

Temperature	25°C	30°C ⁽¹⁾	37°C ⁽¹⁾
Men	6-28 U/l	8-38 U/l	11-50 U/l
Women	4-18 U/l	5-25 U/l	7-32 U/l

⁽¹⁾Calculated

PROCEDURE LIMITATIONS

See Known Interfering Substances under SAMPLE.

PERFORMANCE

a) Reproducibility: CLSI protocol EP15-A was applied. Two activity levels were tested, in replicates by four, during 5 days. With the obtained data, total and intra-assay precision were calculated.

Intra-assay precision (n = 20)

Level	S.D.	C.V.
37.0 U/l	± 0.45 U/l	1.21 %
160.8 U/l	± 0.90 U/l	0.56 %

Total precision (n = 20)

Level	S.D.	C.V.
37.0 U/l	± 1.10 U/l	2.98 %
160.8 U/l	± 3.69 U/l	2.29 %

b) Linearity: manually, the reaction is linear up to 0.200 O.D. $\Delta A/\text{min}$ (250 U/l). For higher values, dilute sample 1/5 or 1/10 with saline solution and repeat assay under the same assay conditions. Multiply results by the dilution performed. In autoanalyzers, a linearity up to 1200 U/l may be observed.

c) Analytical sensitivity: depends on the photometer used and wavelength. In spectrophotometer with 1 cm optical length square cuvettes, for $\Delta A/\text{min}$ of 0.001, the minimum detectable activity change will be of 1 U/l.

PARAMETERS FOR AUTOANALYZERS

For programming instructions check the user manual of the autoanalyzer in use.

WIENER LAB. PROVIDES

100 ml: - 4 x 20 ml Reagent A
- 1 x 20 ml Reagent B
(Cat. 1421404)

100 ml: - 4 x 20 ml Reagent A
- 1 x 20 ml Reagent B
(Cód. 1009330)

100 ml: - 4 x 20 ml Reagent A
- 1 x 20 ml Reagent B
(Cód. 1009258)

125 ml: - 5 x 20 ml Reagent A
- 2 x 12,5 ml Reagent B
(Cód. 1009616)

REFERENCES

- Lum, G.; Gambino, S.R. - Clin. Chem. 18:358 (1972).
- Burrows, S.; Feldman, W.; McBride, F. - Am. J. Clin. Path. 64/3:311 (1975).
- Szasz, G. - Clin. Chem. 15/2:124 (1969).
- S.S.C.C. - Scand. J. Clin. Lab. Invest. 36:119 (1976).
- Szasz G.; Weinmann G.; Staler F.; Wahlefeld W.; Persijn J. - Z. Klin. Chem. Biochem. 12:229 (1974).
- I.F.C.C. - J. Clin Chem. Clin. Biochem. 21: 633 (1983).
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACCC Press, 5th ed., 2000.
- CLSI: Clinical and Laboratory Standards Institute (ex-NCCLS) - Protocols EP-15A, 2001 / EP-17A, 2004.

Symbols

The following symbols are used in packaging for Wiener lab. diagnostic reagent kits.



This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices



Authorized representative in the European Community



"In vitro" diagnostic medical device



Contains sufficient for <n> tests



Use by



Temperature limitation (store at)



Do not freeze



Biological risks



Volume after reconstitution



Contents



Batch code



Manufactured by:



Harmful



Corrosive / Caustic



Irritant



Consult instructions for use



Calibrator



Control



Positive Control



Negative Control



Catalog number

Wiener Laboratorios S.A.I.C.
Riobamba 2944
2000 - Rosario - Argentina
<http://www.wiener-lab.com.ar>
Dir. Téc.; Viviana E. Cétola
Bioquímica
A.N.M.A.T. Registered product
PM-1102-19



Wiener lab.

2000 Rosario - Argentina

UR120831