Apo Calibrator

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Apo

Calibrator

For the calibration of the apolipoprotein A-I and B determination by immunoturbidimetric methods

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Nr kat. 1999722

USES

Apo Calibrator Turbitest AA has been designed for the calibration of apolipoprotein A-I and B determination using Wiener lab's immunoturbidimetric kits **Apo A-I Turbitest AA** and **Apo B Turbitest AA**.

Since it is lot-specific, check the given value on the label.

PROVIDED REAGENTS

Calibrator: lyophilized human serum containing apolipoprotein A-I and B.

NON-PROVIDED REAGENTS

Distilled or demineralized water.

INSTRUCTIONS FOR USE

Reconstitute adding the quantity of distilled or demineralized water stated on the label.

Cap, let stand during 30 minutes and homogenize gently by inversion until complete dissolution before use.

It is recommended to keep reagent inside the original bottle after reconstitution and while is being used.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Avoid ingestion and contact with the skin.

This calibrator has been prepared from non-reactive material to HBsAg, HCV and HIV. However, it should be handled as if it were an infectious material.

Use reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Calibrator: is stable at 2-10°C until the expiration date stated on the box. Once reconstituted it is stable for 1 month at 2-10°C or for 3 months at -20°C aliquoted and frozen immediately after reconstitution. Thawing is recommended only once.

Keep the bottle tightly closed after use. Avoid contamination.

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Wiener Laboratorios S.A. I.C.