SIGNIFICACION CLINICA

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ASO-látex

directo

Slide test for determination of antistreptolysin O

SUMMARY

Since streptococcal infections are frequent, anti-streptolysin O (ASO) antibody is widely found in human in low titers. However, high ASO titer or increasing in ASO titer indicate recent infection produced by group A beta hemolytic Streptococci, such as tonsillitis, scarlet fever, puerperal sepsis, erysipelas, etc.

PRINCIPLE

Antistreptolysin O antibodies are detected in serum when they react to the streptolysin O adsorbed onto an inert latex support. Antistreptolysin antibodies react with the streptolysin producing a macroscopically visible agglutination.

PROVIDED REAGENTS

A. Reagent A: suspension of polystyrene latex particles coated with streptolysin O.

Positive Control*: serum containing antistreptolysin O in concentrations over 250 UI/ml.

Negative Control*: dilution of non-reactive serum proteins.

NON-PROVIDED REAGENTS

Saline solution (for the semi-quantitative technique).

INSTRUCTIONS FOR USE

Reagent A: shake well to homogenize and then replace the screw cap with the provided dropper cap.

Positive Control: ready to use.

Negative Control: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

The Positive and Negative Controls have been tested for the Hepatitis B virus surface antigen, Hepatitis C virus and antibodies against HIV 1/2 and have been found non-reactive. However, they should be handled as infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents are stable in refrigerator (2-10°C) until the expiration date shown on the box. Do not freeze.

INSTABILITY OR DETERIORATION OF REAGENTS

Self-agglutination of the Reagent A is a sign of its deterioration. Discard in that case.

SAMPLE

Serum

a) Collection: obtain serum in the usual way.

b) Additives: not required.

c) Known interfering substances: lipemic or contaminated sera may yield false positive results.

d) Stability and storage instructions: antistreptolysin in serum is stable for 24 hours in refrigerator (2-10°C) or frozen for long periods.

REQUIRED MATERIAL

1- Provided

- plastic or glass slides with black background*
- 25 ul dropper cap

2- Non-provided

- pipettes and micropipettes for measuring the stated volumes
- disposable mixing rods
- stopwatch
- lamp or light source
- Kahn tubes

PROCEDURE

Bring the reagents and the sample to room temperature. Shake the Reagent A before use, after emptying the dropper pipette.

I- QUALITATIVE TECHNIQUE

In one of the circles of the provided slide, place:

Reagent A (25 ul)	1 drop
Sample or Control	25 ul

Mix with a disposable rod (one for each sample) to obtain a uniform suspension all over the circle surface. Immediately, start the stopwatch, gently move the slide and observe the result macroscopically within 2 minutes.

II- SEMI-QUANTITATIVE TECHNIQUE (Titration)

Positive sera can be titrated performing serial dilutions in Kahn tubes.

a) Place 0.5 ml Saline Solution in each tube.

b) Add 0.5 ml serum to tube N° 1 and mix.

Transfer 0.5 ml of this dilution to tube N° 2 and mix; continuing with the dilutions in this manner up to the last tube. The obtained dilutions are equal to 1:2, 1:4, 1:8, 1:16, etc.

c) Test each dilution following the Qualitative Technique.

INTERPRETATION OF RESULTS

Negative: homogeneous suspension.

Positive: agglutination that appears within 2 minutes. It is qualified from 1 to 4 (+).

Titer: it is the inverse of the highest dilution showing a macroscopically visible agglutination.

The approximate level of Antistreptolysin O in the sample, can be calculated by the following formula:

ASO (UI/ml) = Titer x Reaction Sensitivity (200 UI/ml)

Example: the sample has a titer of 1:2. ASO level is:
 $2 \times 200 = 400 \text{ UI/ml}$

QUALITY CONTROL METHOD

Simultaneously process the provided controls, or tested reactive sera, using 25 ul of the corresponding control and 25 ul of Reagent A, following the qualitative technique.

REFERENCE VALUES

Up to 200 UI/ml

It is recommended that each laboratory establishes its own reference values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

- Reaction times longer than 2 minutes may produce false positive reactions as the reagents dry.
- In children older than 6 months and younger than 6 years or with a recent infection, false negative results may be obtained.
- In this kind of infections, an isolated result is only an additional information. Thus, it is recommended to perform serial determinations every 15-20 days for 4 or 6 weeks.
- Any alteration in the Sample/Reagent ratio may produce erroneous results.

PERFORMANCE

a) Sensitivity: 200 UI/ml.

b) Specificity: in an adult healthy population, 95% of the anti-streptolysin O titers obtained are below or equal to 200 UI/ml. In an infant population (over 2 years old) titers up to 300 UI/ml may be obtained.

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Wiener Laboratorios S.A.I.C.