Carbon Dioxide

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Carbon Dioxide

For carbon dioxide (CO₂) determination in serum or plasma

SUMMARY

CO₂ testing in serum represents a measure of serum bicarbonate (HCO₃⁻) in serum. In the body, 95% of CO₂ is present as HCO₃⁻; thus most of what is measured in the laboratory is HCO₃⁻.

CO₂ levels in blood are regulated by the kidneys and lungs (respiratory function).



HCO₃⁻ levels above normal may indicate excessive vomiting, respiratory distress (respiratory conditions), hyperaldosteronism or Cushing syndrome.

HCO₃⁻ levels below normal, may indicate ketoacidosis, lactic acidosis, renal disease, diarrhea, methanol poisoning, toxicity by salicylates (such as aspirin overdose), ethylene glycol poisoning or Addison's disease (adrenal insufficiency), among other pathological conditions.

PRINCIPLE

The method is based on two coupled enzymatic reactions involving phosphoenolpyruvate carboxylase (PEPC) and malate dehydrogenase (MDH).

PEPC catalyzes the first reaction that produces oxaloacetate. In the presence of MDH, NADH is oxidized by oxaloacetate. This results in a decrease in absorbance at 405 or 415 nm which is directly proportional to the concentration of CO₂ in the sample.

PROVIDED REAGENTS

A. Reagent A: PEP 6 mmol/L, PEPC > 100 UI/L, NADH 0.4 mmol/L and MDH > 1000 UI/L in Tris buffer 25 mmol/L. Contains 0.05% sodium azide.

NON-PROVIDED REAGENTS

Carbon Dioxide Calibrator from Wiener lab.

INSTRUCTIONS FOR USE

Reagent A: ready to use.

WARNING

The reagents are for "in vitro" diagnostic use.

Do not pipette by mouth.

Use the reagents according to the working procedures for clinical laboratories.

All reagents and samples must be discarded according to local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date shown on the box

INSTABILITY OR DETERIORATION OF REAGENTS

Reagent A must be a clear solution. If turbid, the reagent may be deteriorated.

Avoid contamination of the reagent with CO₂.

Do not blow into the pipette, because breathing is high in CO₂. Do not leave open containers unnecessarily, since the CO₂ in the air could contaminate the reagent.

Keep container tightly closed.

SAMPLE

Serum or heparinized plasma

a) Collection: obtain the sample in the usual way. Preferably, collect venous blood anaerobically.

b) Additives: in case plasma is used, heparin may be used as anticoagulant for collection. Do not use citrate or oxalate.

c) Known interfering substances: no significant interference have been observed with ascorbic acid up to 5,0 mg/dL, free and conjugated bilirubin up to 40 mg/dL, hemoglobin up to 400 mg/dL and triglycerides up to 1000 mg/dL.

Refer to the bibliography of Young for the effects of drugs in this method.

d) Stability and storage instructions: plasma or serum should be rapidly separated from the clot or cells. Store in tightly closed tubes. CO₂ in blood is stable for 1 hour at 2-10°C anaerobically or for 6 months frozen (-20°C). Centrifuge samples containing precipitate before performing the assay.

REQUIRED MATERIAL (non-provided)

- Automatic analyzer

- Micropipettes and pipettes for measuring the stated volumes.

PARAMETERS FOR AUTOMATIC ANALYZERS

Primary wavelength..... 405 or 415 nm

Reaction type..... kinetic

Reaction direction..... decrease

Reaction temperature..... 37°C

Sample or Calibrator volume..... 3 uL

Reagent A volume..... 300 uL

Incubation time for Reagent A..... 120 sec

Incubation time for Sample..... 120 sec

Reading time..... 60 sec

PROCEDURE

(Automatic analyzers)

Detailed below there is a general procedure for **Carbon Dioxide** in an automatic analyzer. When implementing the technique for a particular analyzer, follow its work instructions.

Reagent A	300 μ L
------------------	-------------

Incubate 120 seconds at 37°C.

Sample or Calibrator	3 μ L
-----------------------------	-----------

Incubate 120 seconds at 37°C. Read kinetically at 405 or 415 nm for 60 seconds.

REFERENCES

- Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982)
- Contarow and Trumper, Clinical Biochemistry, 7th ed., Al Latner, Editor, Saunders, Philadelphia, p. 399 (1975)
- Clinical Chemistry, LA Kaplan, AJ Pesce, Editors, CV Mosby Company, St. Louis (MO), p. 1056 (1984)
- Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, p. 2181 (1994)
- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed. (2007).

CALIBRATION

Use **Carbon Dioxide Calibrator** from Wiener lab, according to the instructions to calibrate the procedure. Calibrator concentrations may vary batch to batch. Enter the concentration value of the calibrator each time a batch is changed.

QUALITY CONTROL METHOD

Carbon Dioxide Control from Wiener lab.

REFERENCE VALUES

Adults: 22-28 mmol/L

It is recommended that each laboratory establish its own reference values.

PROCEDURE LIMITATIONS

See Known Interfering Substances under SAMPLE
Samples with carbon dioxide levels higher than 50 mmol/L should be diluted with 0.9% saline and retested incorporating the dilution factor to calculate the value.

PERFORMANCE

a) Reproducibility: assessed according to EP5-A protocol from CLSI. This study used two samples tested with 2 daily runs in duplicate for 20 days.

Intra-assay precision (n = 80)

Level	S.D.	C.V.
24.1 mmol/L	± 0.56 mmol/L	2.3%
40.1 mmol/L	± 0.91 mmol/L	2.3%

Inter-assay precision (n = 80)

Level	S.D.	C.V.
24.1 mmol/L	± 0.68 mmol/L	2.8%
40.1 mmol/L	± 1.32 mmol/L	3.3%

b) Detection limit: 1.12 mmol/L.

c) Linearity: the reaction is linear up to 50 mmol/L.

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- 1 x 20 ml (Nr kat. 1009915)

ŹRÓDŁA

- Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982).
- Contarow and Trumper, Clinical Biochemistry, 7th ed., Al Latner, Editor, Saunders, Philadelphia, p. 399 (1975).
- Clinical Chemistry, LA Kaplan, AJ Pesce, Editors, CV Mosby Company, St. Louis (MO), p. 1056 (1984).
- Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, p. 2181 (1994).
- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed. (2007).

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Wiener Laboratorios S.A.I.C.