CMD Cleaning Solution

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Nr kat. 1009661

USE

CMD Cleaning Solution is intended to be used in automated clinical chemistry analyzers (Wiener lab. CMD 800 series) for the cleaning process of the electrodes flowcell and pipeline of the ISE module.

PROVIDED REAGENTS

CMD Cleaning Solution: sodium hypochlorite 3%.

INSTRUCTIONS FOR USE

Ready to use.

WARNINGS

- Reagent is for "in vitro" diagnostic use.
- CMD Cleaning Solution is corrosive. H315+H320: Causes skin and eye irritation. H314: Causes severe skin burns and eye damage. P262: Do not get in eyes, on skin, or on clothing. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P280: Wear protective gloves/protective clothing/eye protection/face protection.
- Use the reagent according to the working procedures for clinical laboratories.
- The product should be discarded according to the local regulations in force

STABILITY AND STORAGE INSTRUCTIONS

CMD Cleaning Solution: stable at room temperature (2-30°C) until de expiration date stated on the box.

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Wiener Laboratorios S.A.I.C.