Coagulation

Control P

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- 2 x 1 ml (Cód. 1937005)

REFERÊNCIAS

- Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline. EP15-A2, Vol. 25, Nº 17, NCCLS.
- User Protocol for Evaluation of Qualitative Test Performance; Approved Guideline. EP12-A, Vol. 22, Nº 14, NCCLS.
- Evaluation of Stability of In Vitro Diagnostic Reagents. Approved Guideline. EP-25A, Vol. 29, Nº 20, CLSI (ex NCCLS).
- Tietz Textbook of Clinical Chemistry, Saunders Co., 3rd ed. (1999).

Coagulation

Control P

For precision control in coagulation

USE

Coagulation Control P is designed for use as control in screening tests (prothrombin time, activated partial thromboplastin time, thrombin time, fibrinogen) and in the determination of several coagulation factors with normal and pathological coagulation activity.

See the table of assigned values, since they are lot-specific.

PROVIDED REAGENTS

Control P: pool of pathologic human plasmas with decreased coagulation factors, obtained with citrate as anticoagulant, lyophilized.

NON-PROVIDED REAGENTS

- Distilled water

INSTRUCTIONS FOR USE

- Carefully open the vial to prevent loss of the lyophilized material.
- Add 1.0 ml of distilled water
- Cap and mix by gentle rotation until complete reconstitution, avoiding foaming.
- Dissolve for 10 minutes at room temperature.
- Mix by inversion immediately before use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

The Control was prepared from non-reactive material for HBsAg, antibodies to HCV and antibodies to HIV 1+2. Nonetheless, the calibrator and all blood samples should be handled as if capable of transmitting infections.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagent: stable in refrigerator (2-10°C) until the expiration date shown on the box.

Once reconstituted the reagent is stable for 4 hours at room temperature, for 8 hours in refrigerator (2-10°C) or for 1 month frozen (-20°C). The reagent should be kept frozen only once. Keep the bottles tightly closed.

PROCEDURE

The reconstituted Control should be used in the same way as an unknown sample, following the instructions of the test in use.

PROCEDURE LIMITATIONS

Failures in reconstitution or preservation may cause erroneous results.

See Procedure Limitations in the corresponding Instruction Manual of the kit in use.

WIENER LAB PROVIDES

- 2 x 1 ml (Cat. N° 1937005)

REFERENCES

- Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline. EP15-A2, Vol. 25, N° 17, NCCLS.
- User Protocol for Evaluation of Qualitative Test Performance; Approved Guideline. EP12-A, Vol. 22, N° 14, NCCLS.
- Evaluation of Stability of In Vitro Diagnostic Reagents. Approved Guideline. EP-25A, Vol. 29, N° 20, CLSI (ex NCCLS).
- Tietz Textbook of Clinical Chemistry, Saunders Co., 3rd ed. (1999).

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Wiener Laboratorios S.A.I.C.