Cystatin C

Control

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Cystatin C

Control

For precision control of cystatin C determination

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Nr kat. 1999797

USES

Cystatin C Control Turbitest AA has been designed to be used in the quality control of Wiener lab's **Cystatin C Turbitest AA** kit.

Check the given value chart for the constituents, since they are lot specific.

PROVIDED REAGENTS

Control level 1: human serum containing recombinant human cystatin C in low levels or near to normal values.

Control level 2: human serum containing recombinant human cystatin C in high levels.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

It is recommended to process control daily basis together with patients' samples and set control ranges to suit the needs of each laboratory.

Results should be located within the defined ranges.

Each laboratory should establish guidelines to take corrective actions if values are located outside the interval.

Each laboratory should establish quality control procedures to comply with existing regulations or with accreditation requirements.

WARNINGS

The reagents are for "in vitro" diagnostic use.

Avoid contact with skin and eyes.

Controls have been prepared from non-reactive material for HBsAg, HCV and HIV. All patient samples and controls should be handled as capable of transmitting infectious diseases.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date stated on the box. Once opened, it is stable for up to 4 weeks store at 2-10°C.

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Wiener Laboratorios S.A.I.C.