D-Dimer

Control

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(Cód. 1705025)

REFERÊNCIAS

- Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline. EP15-A2, Vol. 25, Nº 17, NCCLS.
- User Protocol for Evaluation of Qualitative Test Performance; Approved Guideline. EP12-A, Vol. 22, Nº 14, NCCLS.
- Evaluation of Stability of In Vitro Diagnostic Reagents. Approved Guideline. EP-25A, Vol. 29, Nº 20, CLSI (ex NCCLS).
- Tietz Textbook of Clinical Chemistry, Saunders Co., 3rd ed. (1999).

D-Dimer

Control

For precision control in D-dimer determination

USE

D-Dimer Control is designed for precision control of the **D-Dimer** kit from Wiener lab. Refer to the table of assigned values, since they are lot-specific.

PROVIDED REAGENTS

Control level 1: lyophilized human serum containing normal or slightly pathological level of D-Dimer.

Control level 2: lyophilized human serum containing pathological level of D-Dimer.

NON-PROVIDED REAGENTS

Demineralized or distilled water.

INSTRUCTIONS FOR USE

- Open the vials carefully withdrawing the rubber cap to avoid lyophilized material loss.
- Add the quantity of demineralized water stated on the label.
- Cap and mix gently by inversion for 30 minutes, avoiding foam formation.

WARNINGS

Controls are for "in vitro" diagnostic use. Avoid ingestion and contact with skin and eyes.

Controls have been prepared from non-reactive material for hepatitis B surface antigen (HBsAg), antibodies against hepatitis C virus (HCV) and antibodies against human immunodeficiency virus (HIV). However, the Controls and samples should be handled as if they were infective material. Use the reagents keeping the usual work precautions in the clinical analysis laboratory.

All reagents and samples must be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable in refrigerator (2-10°C) until the expiration date indicated on the box.

Once reconstituted they are stable for 24 hours at 2-10°C or for 4 weeks frozen (-20°C).

INSTABILITY OR DETERIORATION OF REAGENTS

Reagents must be clear. Discard if turbidity is observed.

D-Dimer kit. Remove from refrigerator, reconstitute and mix carefully since vigorous agitation can create bubbles. Transfer aliquots of the required volume to the cuvette of the instrument to be used. Cover and store at 2-10°C.

WIENER LAB. PROVIDES

4 x 1 ml: 2 x → 1 ml level 1

2 x → 1 ml level 2

(Cat. No. 1705025)

REFERENCES

- Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline. EP15-A2, Vol. 25, N° 17, NCCLS.
- User Protocol for Evaluation of Qualitative Test Performance; Approved Guideline. EP12-A, Vol. 22, N° 14, NCCLS.
- Evaluation of Stability of In Vitro Diagnostic Reagents. Approved Guideline. EP-25A, Vol. 29, N° 20, CLSI (ex NCCLS).
- Tietz Textbook of Clinical Chemistry, Saunders Co., 3rd ed. (1999).

PROCEDURE

Controls should be used in the same way as an unknown sample, according to the instructions provided in the

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Wiener Laboratorios S.A.I.C.