# D-Dimer

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(Cód. 1705026)

## BIBLIOGRAFIA

- Armando Tripodi, "D-Dimer testing in Laboratory Practice" Clinical Chemistry 57:9, 1256-1262 (2011).
- B. H. Mavromatis, C.M Kessler, "D-Dimer testing: the role of the clinical laboratory in the diagnosis of pulmonary embolism" J. Clin. Pathol. 2001; 54:664-668.
- W.A.M. Lucassen, "Qualitative point of care D-Dimer testing compared with quantitative D-dimer testing in excluding pulmonary embolism in primary care" Journal of Thrombosis and Haemostasis 13: 1004-1009.
- M. Di Nisio, "Diagnostic accuracy of D-Dimer test for exclusion of venous thromboembolism: a systematic review" Journal of Thrombosis and Haemostasis 5; 296-304 (2006).
- Tripodi Armando, "Performance of quantitative D-Dimer methods: results of the Italian external quality assessment scheme". Journal of Thrombosis and Haemostasis 5, 184-185.
- Curtin N., Highe G. Harris, "Extensive evaluation of the IL test D-Dimer immunoturbidimetric assay on the ACL 9000. determines de D Dimer cut off value for reliable exclusion of venous thromboembolism" Lab. Hematol. 2004, 10(2):88-94.

- John D. Olson, "D-Dimer: simple test, tough problems". Archives of Pathology & Laboratory Medicine Aug 2013, Vol. 137, No. 8 (August 2013) pp. 1030-1038.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.
- EP5-A2 (Vol. 24 - N° 25) Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition - CLSI.
- EP6-A (Vol. 23 - N° 16) Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline CLSI.
- EP17-A (Vol.24 - N°34) Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline - CLSI.

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(Cód. 1705026)

## REFERÊNCIAS

- Armando Tripodi, "D-Dimer testing in Laboratory Practice" Clinical Chemistry 57:9, 1256-1262 (2011)
- B. H. Mavromatis, C.M Kessler, "D-Dimer testing: the role of the clinical laboratory in the diagnosis of pulmonary embolism" J. Clin. Pathol. 2001;54:664-668.
- W.A.M. Lucassen, "Qualitative point of care D-Dimer testing compared with quantitative D-dimer testing in excluding pulmonary embolism in primary care" Journal of Thrombosis and Haemostasis 13:1004-1009.
- M. Di Nisio, "Diagnostic accuracy of D-Dimer test for exclusion of venous thromboembolism: a systematic review" Journal of Thrombosis and Haemostasis 5; 296-304 (2006).
- Tripodi Armando, "Performance of quantitative D-Dimer methods: results of the Italian external quality assessment scheme". Journal of Thrombosis and Haemostasis 5, 184-185.
- Curtin N., Highe G. Harris, "Extensive evaluation of the IL test D-Dimer immunoturbidimetric assay on the ACL 9000 determines de D Dimer cut off value for reliable

exclusion of venous thromboembolism" Lab. Hematol. 2004, 10(2):88-94.

- John D. Olson, "D-Dimer: simple test, tough problems". Archives of Pathology & Laboratory Medicine Aug 2013, Vol. 137, No. 8 (August 2013) pp. 1030-1038.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACCC Press, 5th ed., 2000.
- EP5-A2 (Vol. 24 - N° 25) Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition - CLSI.
- EP6-A (Vol. 23 - N° 16) Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline CLSI.
- EP17-A (Vol.24 - N°34) Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline - CLSI.

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 Controle Negativo

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# D-Dimer

Immunoturbidimetric method with latex for the quantitative determination of D-dimer

## SUMMARY

The coagulation process culminates in the cleavage of fibrinogen into fibrin monomers, which spontaneously aggregate to fibrin. Then by action of activated factor XIII, the cross-linking of two D domains is produced by covalent bonds resulting in a solid and stable clot. The fibrinolytic system activates the conversion of plasminogen into plasmin which in turn cleaves the fibrin mesh into smaller fragments. These contain a new antigenic determinant given by the cross-linked D domains (D-dimer), resistant to the action of plasmin. The term D-dimer covers a mixture of fragments and complexes of different molecular weight and its increase indicates an exacerbated coagulation.

High levels of D-dimer are associated with thrombotic diseases, as as deep vein thrombosis (DVT), pulmonary thromboembolism (PE), DIC, myocardial infarction, malignant diseases, obstetric complications, pregnancy, surgeries and post-trauma.

In patients in whom a thrombotic disorder is suspected, a negative result, combined with a low clinical probability "pretest" has a high negative predictive value.

## PRINCIPLE

D-dimer present in the sample agglutinates the latex particles coated with anti-D-dimer antibodies. Turbidity caused by the agglutination of the latex particles is proportional to the concentration of D-dimer in the sample and can be measured spectrophotometrically.

## PROVIDED REAGENTS

**A. Reagent A:** HEPES buffer solution 100 mmol/l.

**B. Reagent B:** suspension of latex particles, coated with anti-D-dimer antibodies in HEPES buffer 10 mmol/l.

**C. Reagent C:** phosphate buffer solution 20 mmol/l.

**Calibrator:** lyophilized human plasma added with D-dimer. See the assigned D-dimer value in the label.

## NON-PROVIDED REAGENTS

- Distilled or dionized water
- **D-Dimer Control** from Wiener lab.

## INSTRUCTIONS FOR USE

**Reagents A, B and C:** ready to use.

Prior to use, it is advisable to let the working temperature balance. Reagent B must be homogenized several times by gentle inversion, before use.

**Calibrator:** reconstitute with 1 ml deionized water. Cover and

leave at room temperature 15 - 30 minutes. Before use, homogenize several times by inversion, avoiding foam formation.

## WARNINGS

The reagents are for "in vitro" diagnostic use.

The reagents contain sodium azide. Avoid contact with skin and eyes. In addition, azide can react with copper and lead pipes forming highly explosive metal azides.

All patient samples should be handled as if they were capable of transmitting infection.

The Calibrator has been tested for HBsAg, HCV and HIV, and is not reactive. However, it must be used as if it were infective material.

Use the reagents keeping the usual work precautions in the clinical chemistry laboratory.

All reagents and samples must be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable in refrigerator (2-10°C) until the expiration date indicated on the box. Do not freeze.

Once opened, they are stable for 4 weeks at 2-10°C.

The reconstituted Calibrator is stable for 12 hours at 2-10°C or for 4 weeks at -20°C. Avoid freezing and thawing more than once.

## SAMPLE

Citrated plasma

**a) Collection:** collect blood carefully, avoiding stasis or foam, and place it in a tube with anticoagulant in a 9 + 1 ratio (example: 4.5 ml blood + 0.5 ml anticoagulant). Mix gently. Centrifuge and separate the plasma before 30 minutes. Collection must be performed with plastic syringes.

**b) Additives:** Anticoagulante TP from Wiener lab. can be used to obtain plasma, or 3.8% or 3.2% sodium citrate. Do not use EDTA or heparin.

**c) Known interfering substances:** haemolysed, lipemic or contaminated samples may interfere with the determination. No interference was observed by bilirubin up to 21 mg/dl, triglycerides up to 350 mg/dl, and hemoglobin up to 250 mg/dl. Refer to Young's bibliography for the effects of drugs in the present method.

**d) Stability and storage instructions:** once separated plasma is stable for 8 hours at 15-25°C, for 24 hours at 2-10°C or for 4 weeks at -20°C. If plasma is to be frozen, the time from sample collection should not exceed 4 hours.

## REACTION CONDITIONS

### General parameters for coagulometers

Name of the test	D Dimer
Primary wavelength	405 or 575 nm
Temperature	37°C
Sample volume	25 µl
Reagent A volume	100 µl
Reagent B volume	50 µl
Incubation time	
Reagent A + Sample	30 seconds
Reading time - ΔTme	180 seconds

Request the applications for the analyzers marketed by Wiener lab. Applications not provided by Wiener lab. must be validated.

### CALIBRATION

Dilute the Calibrator in the following proportions: 1/1, 1/2, 1/4, 1/8 and 1/16 using Reagent C as a diluent. Graph the absorbance change registered according to the concentration of each dilution point, applying a "log-log linear" adjustment. D-dimer concentration is determined by interpolation in the calibration curve. In COR Series coagulometers the calibration is performed automatically by the instrument. A new calibration curve must be performed each time a reagent or calibrator batched is changed.

### QUALITY CONTROL METHOD

**D-Dimer Control** from Wiener lab.

The Controls are processed in the same way as the samples.

### REFERENCE VALUES

<0.40 µg/mL EFU

It is recommended that each laboratory establish its own reference value.

### Unit conversion table

ng/ml DDU	µg/ml or mg/l DDU	µg/ml or mg/l FEU
200	0.20	0.40

DDU: D-dimer units

EFU: Equilavent Fibrinogen Units

### PROCEDURE LIMITATIONS

See Known Interfering substances under SAMPLE.

A new calibration must be performed each time a reagent lot is changed, or when the results of the internal quality control are outside the range considered acceptable. It is recommended that the calibration be repeated every 2 months.

Lipemic samples may yield falsely elevated results.

The presence of rheumatoid factor and HAMA antibodies may cause overestimation of D-dimer values.

D-dimer results should be interpreted together with other clinical

data and complemented with some diagnostic image. The distilled water used for Calibrator reconstitution must be free from preservatives and microbial contamination, otherwise unreliable results may be obtained.

### PERFORMANCE

#### a) Precision:

Samples	Media µg/ml EFU	S.D. <sub>ur</sub>	C.V. <sub>ur</sub> (%)	S.D. <sub>T</sub> µg/ml EFU	CV <sub>T</sub> (%)
Level 1	0.807	0.062	7.7	0.089	7.6
Level 2	1.934	0.062	4.6	0.088	4.6

**b) Sensitivity:** the detection limit is 0.16 µg/mL EFU.

**c) Measuring range:** 0.25 - 12.0 µg/mL EFU.

Samples above 12.0 µg/mL EFU should be diluted 1 + 3 with Reagent C and retested. Multiply the result by 4.

**d) Prozone effect:** no prozone effect is evidenced up to 120 µg/ml EFU.

### WIENER LAB PROVIDES

40 determinations:

2 x 2.5 ml Reagent A

2 x 1 ml Reagent B

1 x 6 ml Reagent C

1 x → 1 ml Calibrator

(Cat. No. 1705026)

### REFERENCES

- Armando Tripodi, "D-Dimer testing in Laboratory Practice" *Clinical Chemistry* 57:9, 1256-1262 (2011).
- B. H. Mavromatis, C.M Kessler, "D-Dimer testing: the role of the clinical laboratory in the diagnosis of pulmonary embolism" *J. Clin. Pathol.* 2001; 54:664-668.
- W.A.M. Lucassen, "Qualitative point of care D-Dimer testing compared with quantitative D-dimer testing in excluding pulmonary embolism in primary care" *Journal of Thrombosis and Haemostasis* 13: 1004-1009.
- M. Di Nisio, "Diagnostic accuracy of D-Dimer test for exclusion of venous thromboembolism: a systematic review" *Journal of Thrombosis and Haemostasis* 5; 296-304 (2006).
- Tripodi Armando, "Performance of quantitative D-Dimer methods: results of the Italian external quality assessment scheme". *Journal of Thrombosis and Haemostasis* 5, 184-185.
- Curtin N., Highe G. Harris, "Extensive evaluation of the IL test D-Dimer immunoturbidimetric assay on the ACL 9000. determines de D Dimer cut off value for reliable exclusion of venous thromboembolism" *Lab. Hematol.* 2004, 10(2):88-94.
- John D. Olson, "D-Dimer: simple test, tough problems". *Archives of Pathology & Laboratory Medicine* Aug 2013, Vol. 137, No. 8 (August 2013) pp. 1030-1038.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.
- EP5-A2 (Vol. 24 - Nº 25) Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition - CLSI.

- EP6-A (Vol. 23 - N° 16) Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline CLSI.
- EP17-A (Vol.24 - N°34) Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline - CLSI.

## SYMBOLS

The following symbols are used in the packaging for Wiener lab. diagnostic reagents kits.

 This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

 Authorized representative in the European Community

 "In vitro" diagnostic medical device

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
(Nr. Kat.. 1705026)

### REFERENCJE

- Armando Tripodi , "D-Dimer testing in Laboratory Practice" Clinical Chemistry 57:9, 1256-1262 (2011).
- B. H. Mavromatis, C.M Kessler, "D-Dimer testing: the role of the clinical laboratory in the diagnosis of pulmonary embolism" J. Clin. Pathol. 2001; 54:664-668.
- W.A.M. Lucassen, "Qualitative point of care D-Dimer testing compared with quantitative D-dimer testing in excluding pulmonary embolism in primary care" Journal of Thrombosis and Haemostasis 13: 1004-1009.
- M. Di Nisio, "Diagnostic accuracy of D-Dimer test for exclusion of venous thromboembolism: a systematic review" Journal of Thrombosis and Haemostasis 5; 296-304 (2006).
- Tripodi Armando, "Performance of quantitative D-Dimer methods: results of the Italian external quality assessment scheme". Journal of Thrombosis and Haemostasis 5, 184-185.
- Curtin N., Highe G. Harris, "Extensive evaluation of the IL test D-Dimer immunoturbidimetric assay on the ACL 9000. determines de D Dimer cut off value for reliable exclusion of venous thromboembolism" Lab. Hematol. 2004, 10(2):88-94.
- John D. Olson, "D-Dimer: simple test, tough problems". Archives of Pathology & Laboratory Medicine Aug 2013, Vol. 137, No. 8 (August 2013) pp. 1030-1038.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.
- EP5-A2 (Vol. 24 - N° 25) Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition - CLSI.

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Wiener Laboratorios S.A.I.C.