# Fecuntest

## directo

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# Fecuntest

## directo

Direct agglutination test for pregnancy diagnosis

### SUMMARY

Human Chorionic Gonadotropin (hCG) is a glycoprotein produced by the placental trophoblastic cells. hCG secretion starts in the earliest gestational stage, rising until it reaches a peak about 9 weeks after the onset of the first missed menstrual period. Levels are quite variable from one woman to another.

Since hCG production is an indication of placental growth, it is possible to establish a direct relationship between the appearance of hCG and pregnancy diagnosis.

### PRINCIPLE

The sample is added to a latex reagent coated with monoclonal antibodies anti-hCG. If sample contains hCG, it will produce a macroscopically visible agglutination.

### PROVIDED REAGENTS

**A Reagent A:** latex-polystyrene particle suspension coated with anti-hCG antibodies.

**Positive Control:** hCG's buffered solution with preservatives.

**Negative Control:** free hCG's buffered solution with preservatives.

### NON-PROVIDED REAGENTS

Saline solution.

### INSTRUCTIONS FOR USE

**Reagent A:** mix well before use.

**Positive Control:** ready to use.

**Negative Control:** ready to use.

### WARNINGS

Reagents are for "in vitro" diagnostic use.

The Positive and Negative Controls have been tested for the Hepatitis B virus surface antigen, Hepatitis C virus and antibodies against HIV 1/2 and have been found non-reactive. However, they should be handled as infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

### STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** are stable at 2-10°C until the expiration date stated on the box. Do not freeze.

### INSTABILITY OR DETERIORATION OF REAGENTS

Turbidity and/or color's change of the reagents may indicate contamination.

### SAMPLE

Urine

**a) Collection:** urine must be collected in a clean container free from detergents.

For the qualitative technique use urine specimen collected at any time of the day. If the first morning urine specimen is used a greater sensitivity in the diagnosis can be achieved. For the titration use 24 hours urine stored at 2-10°C while it is being collected.

**b) Additives:** not required. Do not use preservatives.

**c) Known interfering substances:** hematuria, proteinuria and/or bacterial contamination may lead to erroneous results. Urine samples with sedimentation or turbidity should be filtered or centrifuged before testing.

See Young, D.S. in References for effect of drugs and diseases on the present method.

**d) Stability and storage instructions:** use fresh urine. Sample can be stored at 2-10°C until the test is performed. Urine must be processed within 12 hours from collection.

### REQUIRED MATERIAL

#### 1- Provided

- Plastic or glass slides with black background
- Disposable plastic mix-droppers

#### 2- Non-provided

- Stopwatch
- Proper volumetric material
- Light source

### PROCEDURE

Bring reagents and samples to room temperature before use. Use the plastic mix-dropper provided, in vertical position, to place each sample or control.

#### I- QUALITATIVE TECHNIQUE

In one of the delimited areas of the provided plate, place:

**Sample or Controls** 1 drop

**Reagent A** 1 drop

With the closed end of the plastic mix-dropper used to place the sample or controls, mix for 4-5 seconds to obtain a uniform solution on the surface of the area. Immediately start stopwatch, swing gently the plate and macroscopically observe the result within 2 minutes. Use a light source located directly above the plate.

## II- TITRATION

Mix the urine collected during 24 hours (see SAMPLE) and measure its volume. Take an aliquot and process as follows:

- a) Prepare serial dilutions (1/2, 1/4, 1/8, etc.) using Saline Solution as a diluent.
- b) Test each dilution according to the QUALITATIVE TECHNIQUE.

sitivity will yield negative results.

Very dilute urine specimens may not contain representative levels of hCG. Repeat the test using the first morning urine as sample.

The following may cause erroneous results:

- Exchange the cap's reagents.
- Not performing the mixing times according to PROCEDURE.
- Not using the droppers in vertical position.
- Not using a glass thoroughly clean.
- Not bringing reagents and samples to room temperature before testing.

## PERFORMANCE

**Sensitivity:** **Fecuntest directo** detects hCG in urine at a concentration of 0.2 IU/ml.

This sensitivity has been verified against the hCG Third International Standard from W.H.O.

## WIENER LAB. PROVIDES

Kit for 50 tests (Cat. Nr. 1323152).

## REFERENCES

- Derman, R. et al - "Current Status of Immunologic Pregnancy Test" - Int. Gynaecol. Obstet. 17/90 (1979).
- Reiss, A.M. - "Immunologic Cross-Reaction between Luteinizing Hormone and Antisera to Human Chorionic Gonadotropin" - Int. Arch. Allergy 30:561 (1966).
- Wide, L. - "An Immunological Method for the Assay of Human Chorionic Gonadotropin" - Acta Endocrinol. 41:49. Suppl. 70 (1962).
- Batzer, F. - Fertil. Steril. 34:1 (1980).
- Catt, K.; Dufan, M.; Vaitukaitis, J. - J. Clin. Endocrinol. Metab. 40:537 (1975).
- Braunstein, G.; Rasor, J.; Adler, D.; Danzer, H.; Wade, M. - Am. J. Obstet. Gynecol. 126:678 (1976).
- Lenton, E.; Neal, M.; Sulaiman, R. - Fertil. Steril. 37:773 (1982).
- Dawood, M.; Saxeba, B.; Landesman, R. - Ob. Gyn. 126:678 (1976).
- Braunstein, G. et al. - Ann. Inter. Med. 78:39 (1973).
- Young, D.S. - Effects of Drugs in Clinical Laboratory Tests, 3<sup>rd</sup> Ed., AACC Press, Washington DC, (1990).
- Young, D.S. - "Effects of Diseases on Clinical Laboratory Tests", AACC Press, 4<sup>th</sup> ed., 2001.

## INTERPRETATION OF RESULTS

**Positive:** when the agglutination appears within two minutes.

**Negative:** homogeneous suspension. After two minutes, nonspecific agglutinations may appear.

**Titer:** is the inverse of the highest dilution showing a macroscopically visible agglutination.

## CALCULATIONS

**1- IU of hCG/ml = S x T**

where:

S =sensitivity of the reaction = 0.2 IU/ml of hCG

T = titer

**2- IU of hCG/24 hours = S x T x V**

where:

V = total volume (ml) of urine excreted during 24 hours.

## QUALITY CONTROL METHOD

**Positive Control:** positive reactions of **Fecuntest directo** show a lumpy appearance which can be confirmed testing a drop of urine containing higher or equal hCG's quantity to the kit's sensitivity or a Positive Control drop instead of sample, following the described technique (see PROCEDURE).

The strongly positive reactions can be stored in aliquots at -20°C for no more than 60 days, and should be discarded after their use, since their integrity does not last after repeated freezing and thawing.

**Negative Control:** negative reactions of **Fecuntest directo** show homogeneous appearance which can be confirmed processing a drop of saline solution or a drop of Negative Control instead of Sample, following the technique previously explained (see PROCEDURE).

## REFERENCE VALUES

During the first trimester of pregnancy, values between 25-50 IU/ml of hCG are observed, reaching a maximum of 250 IU/ml in normal conditions.

A number of conditions, in addition to pregnancy, including trophoblastic disease and certain non-trophoblastic neoplasms produce high hCG levels. Such diagnosis should be considered in case the clinical evidence is appropriate.

A final clinical diagnosis should not be based on one single test. Physician should perform the diagnosis once the clinical and laboratory tests have been assessed.

## PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

Samples with hCG concentrations below the method's sen-

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Wiener Laboratorios S.A.I.C.