Fibrinogen

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Fibrinogen

Immunoturbidimetric method for quantitative determination of fibrinogen in plasma

SUMMARY

Fibrinogen is a glycoprotein present in plasma and platelet alpha-granules. It is the clotting factor found in highest concentration in plasma.

In the presence of trauma or vascular injury, thrombin cleaves fibrinogen making fibrin monomers which spontaneously polymerize and are stabilized leading to insoluble fibrin mesh.

Low fibrinogen levels can be found in hereditary disorders such as hypofibrinogenemia, afibrinogenemia, dysfibrinogenemia and in other circumstances such as liver disease, disseminated intravascular coagulation, fibrinolytic syndromes, etc.

Increased levels are present in diabetes, inflammatory disease, etc.

In addition, high fibrinogen levels are now recognized as a risk factor for the development of cardiovascular disease.

PRINCIPLE

Fibrinogen reacts to the specific antibody forming insoluble immune complexes. The turbidity caused by these immune complexes is proportional to fibrinogen concentration in the sample and may be spectrophotometrically measured.

PROVIDED REAGENTS

A. Reagent A: phosphate buffer, pH 7.4.

B. Reagent B: polyclonal antibodies anti-human fibrinogen (goat) in phosphate buffer, pH 7.4.

NON-PROVIDED REAGENTS

- Saline solution

- Wiener lab.'s **Fibrinogen Calibrator Turbitest AA**

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

All patient samples should be handled as capable of transmitting infectious diseases.

Use the reagents according to the working procedures for clinical laboratories.

Reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date stated on the box. Do not freeze.

SAMPLE

Plasma

a) Collection: obtain blood samples carefully avoiding stasis or foaming, mix gently in a tube with anticoagulant in 9 + 1 proportion (example: 4.5 ml blood + 0.5 ml anti-coagulant).

Centrifuge and remove plasma before 30 minutes. Withdraw the specimen using plastic syringes.

b) Additives: Wiener lab.'s **Anticoagulante TP** or sodium citrate 3.8% or 3.2% could be used to obtain plasma.

c) Known interfering substances: do not use hemolyzed, lipemic or contaminated samples. Before testing, particles in samples should be removed by centrifugation.

No interferences have been observed with hemoglobin up to 1000 mg/dl, triglycerides up to 1700 mg/dl, bilirubin up to 30 mg/dl and rheumatoid factor up to 500 IU/ml.

See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: sample should be preferably fresh. In case it cannot be processed immediately, the sample can be kept for up to 48 hours at 2-10°C or for longer period store at -20°C.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer.

- Square spectrophotometric cuvettes.

- Micropipettes and pipettes for measuring the stated volumes.

- Kahn or hemolysis tubes.

- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 340 nm

- Reaction temperature: room temperature (25°C). Temperature control is not critical, it can range between 22 and 30°C.

- Reaction time: 15 minutes

- Sample volume: 6 ul

- Final reaction volume: 1806 ul

Sample and reagents volumes may be proportionally changed without affecting the calculation factors.

PROCEDURE

CALIBRATION CURVE

In Kahn tubes dilute the **Fibrinogen Calibrator Turbitest AA** with saline solution 1:1, 1:2, 1:4, 1:8 and 1:16, using saline solution as the zero point.

Diluted Calibrator

6 ul

Reagent A	1500 ul
Homogenize and measure the absorbance of each dilution at 340 nm (OD ₁), setting the instrument to zero with distilled water. Then, add:	
Reagent B	300 ul
Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD ₂), setting the instrument to zero with distilled water. Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each Fibrinogen Calibrator dilution, including the zero point. Draw on graph paper the ΔA absorbance differences based on the Fibrinogen Calibrator concentration in mg/dl (g/l).	
SAMPLES PROCEDURE	
Sample	6 ul
Reagent A	1500 ul
Homogenize and measure the absorbance at 340 nm (OD ₁), setting the instrument to zero with distilled water. Then add:	
Reagent B	300 ul
Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD ₂), setting the instrument to zero with distilled water.	

Intra-assay precision

Level	S.D.	C.V.
63.9 mg/dl	± 1.4 mg/dl	2.3%
254.6 mg/dl	± 9.1 mg/dl	3.6%
365.4 mg/dl	± 17.3 mg/dl	4.7%

b) Detection limit: 13 mg/dl

c) Measuring range: 13 - 500 mg/dl

d) Prozone effect: not noted until 5200 mg/dl Fibrinogen.

WIENER LAB. PROVIDES

60 ml: - 1 x 50 ml Reagent A

- 1 x 10 ml Reagent B

(Cat. N° 1009363)

60 ml: - 1 x 50 ml Reagent A

- 1 x 10 ml Reagent B

(Cat. N° 1009957)*

REFERENCES

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.

- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.

CALCULATIONS

Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each sample tested. Interpolate this ΔA in the calibration curve to determine the concentration in mg/dl (g/l) corresponding to the sample under study. Samples with an absorbance above that of the Fibrinogen Calibrator Turbitest AA must be diluted with saline solution and processed again. Multiply the obtained result by the dilution factor.

QUALITY CONTROL METHOD

It is recommended the use of **Plasma Control normal-patológico**, separately provided by Wiener lab.

The Control is processed in the same manner as samples.

REFERENCE VALUES

200 - 400 mg/dl (2 - 4 g/l)

Each laboratory should set its own reference values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

It is recommended to perform a complete recalibration when changing reagent lot or when suggested by Quality Control. Avoid contamination to preserve the integrity of the reagents. Only use thoroughly clean and dry micropipettes for measurement.

PERFORMANCE

a) Reproducibility: replicates of samples containing different fibrinogen levels were assayed and the following results were obtained:

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Wiener Laboratorios S.A.I.C.