HbA1c

enzymatic - Calibrator

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HbA1c

enzymatic - Calibrator

For the calibration of HbA1c determination

USES

HbA1c enzymatic Calibrator is designed for the calibration of HbA1c determination by using **HbA1c enzymatic** kit.

PROVIDED REAGENTS

Calibrator 1: lyophilized human hemoglobin with HbA1c level within the range of normal values or within the clinically expected values for controlled diabetic patients (4-7%)
Calibrator 2: lyophilized human hemoglobin with HbA1c level within the range of pathological values (8-14%).

NON-PROVIDED REAGENTS

Demineralized water

INSTRUCTIONS FOR USE

Open the vial by slowly removing the rubber cap to avoid loss of lyophilized material.

Add the demineralized water quantity indicated on the vial label. Cap and mix by gentle inversion several times for a period of 30 minutes, avoiding foam formation.

Calibrators require pretreatment using **HbA1c enzymatic Lysis Buffer**.

WARNING

Reagents are for "in vitro" diagnostic use.

Do not ingest. Avoid contact with skin and eyes.

Calibrators were prepared from non-reactive material for HBsAg, HCV and HIV. However, calibrators should be handled as if they were infective samples.

Use the reagents according to the working procedures for clinical laboratories.

All reagents and samples should be discarded according to current regulations.

STABILITY AND STORAGE INSTRUCTIONS

Lyophilized calibrators are stable at 2-10°C until the expiration date stated on the box. After reconstitution they are stable for up to 14 days at 2-10°C. After use, keep bottles tightly capped and stored at 2-10°C.

PROCEDURE

Calibrators must be used in the same way as an unknown sample, according to the instructions for use that comes with **HbA1c enzymatic** kit. Remove from refrigerator, reconstitute and mix carefully because vigorous agitation can create bubbles. Transfer aliquots of the required volume to the cuvette of the instrument to be used. Cap and store at 2-10°C.

%HbA1c values may be expressed according to IFCC or DCCT/NGSP standardization using the following calculation methods:

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Wiener Laboratorios S.A.I.C.