Lipid

Control

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- 2 x 3 ml (Cód. 1999740).

REFERÊNCIAS

- Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982)
- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed. (2007).

Lipid

Control

For precision control of lipid determination

USES

Lipid Control is a liquid control designed to control the precision of lipid determination in serum or plasma. Check the given value chart for the two levels provided, since they are lot specific.

PROVIDED REAGENTS

Control (1 and 2): liquid controls prepared from human serum containing two different levels of lipids.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

Homogenized the content of each vial by inversion before testing. Do not shake.

Remove the amount of control to use. Do not return the used material into the original vial.

After use, cap immediately and store refrigerated (2-10°C). Avoid controls contamination.

WARNINGS

Reagents are for "in vitro" diagnostic use.

These control sera have been prepared from non-reactive material to HBsAg, HCV and HIV. However, the controls and samples should be handled as if it were infectious material. Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable in refrigerator at 2-10°C until the expiration date stated on the box.

INSTABILITY OR DETERIORATION OF REAGENTS

Any variation of the organoleptic features or contamination of the controls may be a sign of its deterioration. The presence of turbidity or microbial growth is evidence of contamination and deterioration of control, so it should be discarded.

PROCEDURE

Controls should be used in the same manner as an unknown sample, following the package insert of the reagent kit in use.

REFERENCE VALUES

The obtained values should be within the acceptable range

established. However it is recommended that each laboratory establish its own mean and use the provide only as a guideline.

PROCEDURE LIMITATIONS

See PROCEDURE LIMITATIONS on the package insert of the kit in use.

The use of different methods for Wiener lab kits. to establish the expected values can lead to different values than stated.

WIENER LAB. PROVIDES

- 2 x 3 ml (Cat. 1999740).

REFERENCES

- Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982)
- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed. (2007).

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Wiener Laboratorios S.A.I.C.