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PCR-látex

directo

Direct latex agglutination test for the determination of C Reactive Protein

SUMMARY

C Reactive Protein (CRP) is a thermolabile protein that does not cross the placental barrier and has an electrophoretic mobility between the α and β globulin's zones. Its name comes from its ability to precipitate the pneumococci C-polysaccharides. It is one of the acute phase proteins, and it increases in serum in a large variety of inflammatory diseases or in response to tissue necrosis.

Its determination is important since it rapidly increases at the onset of the disease, 14 to 26 hours after the tissue inflammation or injury, and disappears during the recovery stage. It only appears during the active phase of the inflammatory process. CRP is usually increased in active rheumatoid arthritis, viral infections, tuberculosis, active rheumatic fever, acute myocardial infarction, etc. It can also be found after surgery and in a large proportion after blood transfusions. CRP determination not only indicates the intensity of the disease, but also the patient's response to a given therapy.

PRINCIPLE

CRP is detected in serum by reacting to a specific antibody adsorbed onto inert latex particles. CRP binds to the adsorbed antibodies producing latex particles agglutination.

PROVIDED REAGENTS

A. Reagent A: suspension of latex-polystyrene particles sensitized with anti-CRP antibodies.

Positive Control*: serum protein dilution containing C reactive protein.

Negative Control*: non-reactive serum protein dilution.

NON-PROVIDED REAGENTS

Saline solution.

INSTRUCTIONS FOR USE

Reagent A: mix thoroughly and then change the cap for the additionally provided dropper.

Negative Control: ready to use.

Positive Control: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

The Controls have been tested for the Hepatitis B virus surface antigen, antibodies to Hepatitis C virus and HIV 1/2 and have been found non-reactive. However, they should be handled as infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided reagents are stable in refrigerator (2-10°C) until the expiration date shown on the box. Do not freeze.

SAMPLE

Serum

a) Collection: obtain serum in the usual way.

b) Additives: not required.

c) Known interfering substances: distinctly lipemic or contaminated sera may yield false positive results.

d) Stability and storage instructions: serum should be preferably fresh. If it cannot be processed immediately, it can be kept refrigerated (2-10°C) up to 24 hours and frozen (-20°C) up to 4 weeks.

REQUIRED MATERIAL

1- Provided

- plastic or glass slides with black background*
- one 25 ul dropper cap

2- Non-provided

- Pipettes and micropipettes for measuring indicated volumes.
- Disposable mixing rods.
- Stopwatch.
- Lamp or light source.

PROCEDURE

Bring the reagents and the samples to room temperature before use. Shake the Reagent A before using, after emptying the dropper pipette.

I- QUALITATIVE TECHNIQUE

Sample or Control	25 ul
Reagent A	1 drop (25 ul)

Mix with provided disposable rod to obtain a uniform suspension all over the circle surface. Immediately start the stopwatch. Gently balance the slide and observe the result macroscopically under a light source, within 2 minutes.

II- TITRATION

Positive sera may be titrated performing serial dilutions in 8 Kahn tubes.

a) Place 0.5 ml saline solution in each tube.

* Non-provided with all kit sizes

b) Add 0.5 ml serum to tube N° 1 and mix.
Transfer 0.5 ml of this dilution to tube N° 2 and mix,
continuing the dilutions up to the last tube. The obtained
dilutions are 1:2, 1:4, 1:8, 1:16, etc.
c) Test each dilution following TECHNIQUE I.

INTERPRETATION OF RESULTS

Negative: homogeneous suspension.

Positive: agglutination that appears within 2 minutes. It is classified from 1 to 4 +.

Titer: is the inverse of the highest dilution showing a macroscopically visible agglutination.

The approximate CRP concentration in the sample may be calculated using the following formula:

$CRP (mg/l) = \text{Titer} \times \text{Reaction sensitivity} (6 \text{ mg/l})$

E.g.: the sample has a 1:2 titer. Its CRP concentration is:
 $2 \times 6 = 12 \text{ mg/l}$

QUALITY CONTROL METHOD

Simultaneously process the provided controls or reactive sera using 25 ul of the corresponding Control and 25 ul Reagent A, following the qualitative technique.

REFERENCE VALUES

Up to 6 mg/l.

It is recommended that each laboratory establishes its own reference values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

- Reaction times longer than 2 minutes may produce falsely positive results, due to the drying of the reagents.
- Any alteration in the Sample/Reagent proportion may produce erroneous results.

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Wiener Laboratorios S.A.I.C.