PCR-látex

directo

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PCR-látex

directo

Direct latex agglutination test for the determination of
C Reactive Protein

SUMMARY

C Reactive Protein (CRP) is a thermolabile protein that does not cross the placental barrier and has an electrophoretic mobility between the α and β globulins zones. Its name is due to its ability to precipitate the pneumococci C-polysaccharides. It is one of the "acute phase" proteins, increasing its concentration in serum in a large variety of inflammatory diseases or in response to tissue necrosis.

Its determination is important since it rapidly increases at the onset of the disease, 14 to 26 hours after the tissue inflammation or injury, and disappears during the recovery stage. It only appears during the active phase of the inflammatory process. CRP is usually increased in: active rheumatoid arthritis, viral infections, tuberculosis, active rheumatic fever, acute myocardial infarction, etc. It can also be found after surgery and in a large proportion after blood transfusions.

CRP determination not only indicates the intensity of the disease, but also the patient's response to a given therapy.

PRINCIPLE

CRP is detected in serum by reacting to a specific antibody adsorbed onto inert latex particles. CRP binds to the adsorbed antibodies producing the latex particles agglutination.

PROVIDED REAGENTS

A. Reagent A: suspension of latex-polystyrene particles sensitized with anti-CRP antibodies.

Negative Control: negative serum dilution.

Positive Control: positive serum dilution.

NON-PROVIDED REAGENTS

Saline solution.

INSTRUCTIONS FOR USE

Reagent A: mix thoroughly and then replace screw cap with provided dropper.

Negative and Positive Control: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

The Controls have been tested for the Hepatitis B virus surface antigen, antibodies to Hepatitis C virus HIV 1/2 and have been found non-reactive. However, they should be handled as infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided reagents are stable at 2-10°C until the expiration date shown on the box. Do not freeze.

INSTABILITY OR DETERIORATION OF REAGENTS

Self agglutination of the Reagent A is a sign of its deterioration. In that case, discard.

SAMPLE

Serum

a) Collection: obtain serum in the usual way.

b) Additives: not required.

c) Known interfering substances: distinctly lipemic or contaminated sera may yield false positive results.

d) Stability and storage instructions: serum should be preferably fresh. If it cannot be processed immediately, it can be kept at 2-10°C for up to 24 hours and at -20°C for up to 4 weeks.

REQUIRED MATERIAL

1- Provided

- plastic or glass slides with black background.

2- Non-provided

- Volumetric material for measuring and diluting the samples.

- disposable mixing rods

- Stopwatch.

- Lamp or light source.

PROCEDURE

Bring the reagents and the samples to room temperature before using. Shake the Reagent A before use, previously draining off the dropper pipette.

I- QUALITATIVE TECHNIQUE

Sample	1 drop (50 ul)
Reagent A	1 drop (50 ul)

Mix with a disposable mixing rod to obtain a uniform suspension all over the circle surface, and immediately start the stopwatch. Gently move the slide and observed macroscopically under a source of light, within 2 minutes.

II- TITRATION

Positive sera may be titrated performing serial dilutions in 8 Kahn tubes.

- a) Place 0.5 ml saline solution in each tube.
- b) Add 0.5 ml serum to tube N° 1 and mix. Transfer 0.5 ml of this dilution to tube N° 2 and mix, continuing the dilutions up to the last tube. The obtained dilutions are 1:2, 1:4, 1:8, 1:16, etc.
- c) Test each dilution following I.

INTERPRETATION OF RESULTS

Negative: homogeneous suspension.

Positive: agglutination that appears within 2 minutes. It is qualified from 1 to 4 +.

Titer: is the inverse of the highest dilution showing a macroscopically visible agglutination.

The approximate CRP concentration of the sample can be calculated with the following formula:

CRP (mg/l) = Titer x Reaction sensitivity (6 mg/l)

E.g.: the sample has a 1:2 titer. Its CRP concentration is:
 $2 \times 6 = 12 \text{ mg/l}$

QUALITY CONTROL METHOD

Simultaneously process the provided Negative and Positive controls, using a drop of the corresponding Control instead of the sample.

REFERENCE VALUES

Up to 6 mg/l.

It is recommended that each laboratory establishes its own reference values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

Reaction times longer than 2 minutes may produce false positive results, due to the drying of the reagents.

PERFORMANCE

Sensitivity: PCR-látex *directo* detects 6 mg/l of C Reactive Protein.

WIENER LAB. PROVIDES

Kit for 50 tests (Cat. 1683152).

REFERENCES

- Singer, J.M.; Plotz, C.M.; Parker, E. and Elster, S.K. - Am. J. Clin. Path. 28:611, 1957.
- Nilson, L.A. - Acta Pathol. Microbiol. Scand. 73:129, 1968.
- Scherffarth, F.; Pérez-Miranda, M. and Goetz, H. - Blut. 20:296, 1970.

SYMBOLS

The following symbols are used in the packaging for Wiener lab. diagnostic reagents kits.

This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

	Authorized representative in the European Community
	"In vitro" diagnostic medical device
	Contains sufficient for <n> tests
	Use by
	Temperature limitation (store at)
	Do not freeze
	Biological risks
	Volume after reconstitution
	Contents
	Batch code
	Manufactured by:
	Harmful
	Corrosive / Caustic
	Irritant
	Consult instructions for use
	Calibrator
	Control
	Positive Control
	Negative Control
	Catalog number

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Wiener Laboratorios S.A.I.C.