# Plasma Control

*Normal*

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- 6 x 1 ml (Cód. 1937001).

## REFERÊNCIAS

- Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982)
- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed. (2007).

# Plasma Control

*Normal*

For precision control in coagulation

## USES

The **Plasma Control normal** is designed to be used as control of the determination of prothrombin time, thrombin time, activated partial thromboplastin time and fibrinogen.

## PROVIDED REAGENTS

**Normal Control:** lyophilized pool of normal human plasmas obtained with citrate as anticoagulant.

## NON-PROVIDED REAGENTS

Bidistilled or deionized water.

## INSTRUCTIONS FOR USE

Reconstitute each vial with 1.0 ml bidistilled or deionized water. Allow this material to dissolve for 20 minutes at room temperature, mixing without inversion. Make sure that all particles are completely dissolved.

## WARNINGS

Reagent is for "in vitro" diagnostic use.

The normal control has been prepared from material non-reactive for HBsAg, HCV and HIV. However, the control and all other blood samples should be handled as potentially infectious biological material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

The reagent is stable in refrigerator (2-10°C) until the expiration date indicated on the box. Once reconstituted, it is stable for 8 hours both at room temperature and refrigerated (2-10°C).

## INSTABILITY OR DETERIORATION OF REAGENTS

Any variation in the color of the normal control may be a sign of deterioration.

## PROCEDURE

Reconstituted normal control is used in the same way as an unknown sample, according to the instructions attached to the kit of reagents for prothrombin time, thrombin time, activated partial thromboplastin time and fibrinogen used.

## PROCEDURE LIMITATIONS

Mistakes when reconstituting the reagent may cause erroneous results.

See LIMITATIONS OF THE PROCEDURE on the package insert of the kit in use.

## REFERENCE VALUES

The values for prothrombin time, thrombin time, activated partial thromboplastin time and fibrinogen, using Plasma Control, depend on the reagents and coagulation analyzers in use.

## WIENER LAB. PROVIDES

- 6 x 1 ml (Cat. N° 1937001)

## REFERENCES

- Tietz, N. W. (Ed): Fundamentals of Clinical Chemistry, W. B. Saunders Co., Philadelphia, 865 (1982)
- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed. (2007).

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Wiener Laboratorios S.A.I.C.