Standatrol S-E

2 niveles

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- 6 frascos x 5 ml (Cód. 1937553).

REFERÊNCIAS

- International Federation of Clinical Chemistry. - Clin. Chem. 23/9:1784 (1977).
- Tonks, D. B. - Can. J. Med. Tech. 30:38 (1969).

Standatrol S-E

2 niveles

Lyophilized serum for precision control in clinical chemistry

USES

Standatrol S-E 2 niveles is adaptable to different uses:

- intralaboratory precision studies in which a stable, reproducible in time control specimen is required.
- It can be used in the reproducibility daily control and in the preparation of "control charts" following usual methods;
- as "unknown sample" in interlaboratory precision surveys, based on its constant dosage and high stability.

PRINCIPLE

Standatrol S-E 2 niveles contains the compounds usually determined in clinical chemistry laboratories. It must be noted that the assigned values for the different constituents of the control sera have been obtained with Wiener lab.'s methods and reagents, thus, the obtained results will only be comparable to the list if the corresponding methods and reagents have been used.

PROVIDED REAGENTS

Control Level 1: 5 ml-vials of homogenized and lyophilized serum with normal concentrations or normal/pathological threshold of metabolites and enzymes.

Control Level 2: 5 ml-vials of homogenized and lyophilized serum with pathological concentrations or normal/pathological threshold of metabolites and enzymes.

NON-PROVIDED REAGENTS

Bidistilled or deionized water.

INSTRUCTIONS FOR USE

- Open the vial, slowly pulling out the rubber stopper to avoid any loss of the lyophilized material.
- Add exactly 5.00 ml bidistilled or deionized water (using appropriate volumetric device).
- Cap and mix by gentle inversion, avoiding the formation of foam. Do not shake.
- Let dissolve for 20 minutes at room temperature, occasionally mixing by inversion.
- Mix by inversion immediately before use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

These control sera have been prepared from non-reactive material to HBsAg, HCV and HIV. However, the controls and blood samples should be handled as if it were infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Standatrol S-E 2 niveles: is stable in refrigerator (2-10°C) until the expiration date shown on the box.

Once reconstituted the controls, most of their constituents are stable for 10 days refrigerated (2-10°C), with the exception of alkaline phosphatase whose activity may increase with time, and bilirubin that is stable for 12 hours at 4°C and protected from light. For each particular component, the stability is similar to that in the native serum and may be consulted in the corresponding package insert.

To improve its performance, it is recommended to freeze (-20°C) the reconstituted control serum, divided in fractions in airtight containers (e.g. 1.5 ml microcentrifuge tubes). In this way the controls are stable up to 35 days. They should only be thawed one time at room temperature, homogenizing the serum before use.

INSTABILITY OR DETERIORATION OF REAGENTS

Any variation of the organoleptic features of the control sera (such as unusual coloration, humid appearance of lyophilized reagent, incomplete dissolution with abundant flocculation or lump formation) may be a sign of its deterioration.

PROCEDURE

Reconstituted controls should be used in the same manner as an unknown sample, following the package insert of the reagent kit in use.

Note: the enzymes assigned values were obtained using techniques at 37°C in the different automated analyzers.

PROCEDURE LIMITATIONS

Mistakes in the reconstitution or conservation may cause erroneous results.

See PROCEDURE LIMITATIONS on the package insert of the kit in use.

WIENER LAB. PROVIDES

6 x 5 ml vials (Cat. 1937553).

REFERENCES

- International Federation of Clinical Chemistry - Clin. Chem. 23/9:1784, 1977.
- Tonks, D.B. - Can. J. Med. Tech. 30:38, 1969.

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Wiener Laboratorios S.A.I.C.