



Acid Solution

Acid cleaning solution for use in automated analyzers

USES

Acid Solution is employed in automated analyzers as cleaning solution for cuvettes reading, probes and hydraulic system.

PROVIDED REAGENTS

Acid Solution: hydrochloric acid solution, 1 N.

INSTRUCTIONS FOR USE

Acid Solution: ready to use.

The solution could be used diluted according to the specifications given for each analyzer.

WARNINGS

The reagent is for "in vitro" diagnostic use.

The solution is caustic. R36/38: irritates eyes and skin. S26/28: in case of contact with the eyes and skin, rinse immediately with plenty of water and seek medical advice. S37: wear suitable gloves.

Use the reagents according to the working procedures for clinical laboratories.

All reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Acid Solution: stable at room temperature (< 25°C) until the expiration date stated on the box.

INSTABILITY AND DETERIORATION OF REAGENTS

Evident signs of turbidity or presence of particles in suspension indicate reagent deterioration. Discard in such case. Avoid external contamination.

WIENER LAB. PROVIDES

- 1 x 250 ml (Cat. 1999720).

SYMBOLS


The following symbols are used in packaging for Wiener lab. diagnostic reagents kits.



This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

 Authorized representative in the European Community

 "In vitro" diagnostic medical device

 Contains sufficient for <n> tests

 Use by

 Temperature limitation (store at)

 Do not freeze

 Biological risks

 Volume after reconstitution


 Contents

 Batch code


 Manufactured by:

 Harmful

 Corrosive / Caustic

 Irritant

 Consult instructions for use


 Calibrator

 Control

 Positive Control

 Negative Control

 Catalog number

 Wiener Laboratorios S.A.I.C.
Riobamba 2944
2000 - Rosario - Argentina
<http://www.wiener-lab.com.ar>
Dir. Téc.: Viviana E. Cétola
Bioquímica
A.N.M.A.T. Registered product
PM-1102-44



Wiener lab.

2000 Rosario - Argentina