Ammonia

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- 1 x 20 ml (Cód. 1009911)*

REFERÊNCIA

- Tietz, N.W., Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: W.B. Saunders Co; 1995:44.

- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed., 2007.

- Mondzac A, Ehrlich GE, Seegmiller JE - An enzymatic determination of ammonia in biological fluids - J. Lab. Clin. Med. 66/3:526-31, 1965.

Ammonia


For ammonia determination (NH₃) in plasma

SUMMARY

It has been shown that increased levels of ammonia (hyperammonemia) produce toxic effects on the central nervous system. There are several diseases, both inherited and acquired, causing hyperammonemia. Inherited deficiency of the enzyme of the urea cycle is the leading cause of hyperammonemia in children. Among the acquired causes that produce an increase in ammonia levels are advanced liver disease and kidney failure.

PRINCIPLE

The method is based on the reaction of ammonia in the sample with α -ketoglutarate and NADPH in the presence of glutamate dehydrogenase (GLDH) to yield glutamate and NADP⁺. The decrease in absorbance at 340 nm is proportional to the ammonia concentration in plasma.

PROVIDED REAGENTS

A. Reagent A: 8.0 mmol/L α -ketoglutarate, 0.35 mmol/L NADPH and 1000 kU/L GLDH in 100 mmol/L Tris buffer, pH 8.6.

NON-PROVIDED REAGENTS

Ammonia Calibrator from Wiener lab.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNING

The reagents are for "in vitro" diagnostic use. Use the reagents according to the working procedures for clinical laboratories. All reagents and samples must be discarded according to local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable in refrigerator (2-10°C) until the expiration date shown in the box.

INSTABILITY OR DETERIORATION OF REAGENTS

The reagent solution should be clear. If turbid, the reagent may be deteriorated.

SAMPLE

Plasma with EDTA
(See Procedure Limitations)

a) Collection: collect the sample avoiding stasis and he-

molysis. Centrifuge in a capped tube as soon as possible.

b) Additives: use only EDTA (**Anticoagulante W** from Wiener lab) as anticoagulant.

c) Known interfering substances: visible hemolysis interferes with the assay. No significant interference was observed by bilirubin up to 60 mg/dL, hemoglobin up to 50 mg/dL or triglycerides up to 500 mg/dL. Refer to Young bibliography for the effects of drugs in this method.

d) Stability and storage instructions: place the sample in ice and test immediately. The plasma is stable for 3 hours at 4°C in a closed container. Centrifuge the samples containing precipitate before performing the assay.

REQUIRED MATERIAL (non-provided)

- Micropipettes and pipettes for measuring the stated volumes.
- Spectrophotometer.
- Water bath at the selected reaction temperature.
- Stopwatch.

ASSAY CONDITIONS

Wavelength: 340 nm
Reaction temperature: 25, 30 or 37°C
Reaction time: 5 minutes
Sample volume: 100 μ L
Final reaction volume: 1.1 mL

PROCEDURE

In two cuvettes maintained at the selected temperature, labeled B (blank) and U (sample or calibrator), place:

	B	U
Sample or Calibrator	-	100 μ L
Distilled water	100 μ L	-
Reagent A	1.0 mL	1.0 mL

Mix and pre-incubate for 30 seconds. Read OD₁ absorbance immediately and measure OD₂ after exactly 2.5 minutes.

CALCULATIONS

$$\text{Ammonia (ug/dL)} = \Delta A_M \times \text{factor} \quad \Delta A = \text{OD}_2 - \text{OD}_1$$

$$\text{Factor} = \frac{x \text{ ug/dL}}{\Delta A_C}$$

where:

ΔA_M : Δ Absorbance of the sample

x ug/dL: ammonia concentration in the Calibrator
 ΔA_c : Δ Absorbance of the Calibrator

QUALITY CONTROL METHOD

Ammonia Control from Wiener lab.

REFERENCE VALUES

Men: 25 - 94 ug/dL (14,7 - 55,3 umol/L)

Women: 19-82 ug/dL (11,2 - 48,2 umol/L)

Each laboratory should establish its own reference values. The results should be assessed together with the patient's medical background, medical examination and other laboratory findings.

UNITS CONVERSION TO SI SYSTEM

Ammonia (ug/dL) = 1,703 x ammonia (umol/L)

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

Do not use serum as a sample since during coagulation it can be originate ammonia, producing erroneous results.

Manually dilute samples that are above the measuring range with ammonia-free freshly distilled or deionized water (e.g. 1 +1). Multiply the result by the dilution factor (e.g. 2).

Avoid reporting the results generated from a new automated run unless a fresh sample is added.

PERFORMANCE

a) Reproducibility: evaluated according to EP5-A protocol from CLSI. This study used two different samples of ammonia (NH_3) that were tested performing 2 daily runs by duplicate during 5 days.

Intra-assay precision

Level	S.D.	C.V.
71,9 ug/dL	3,0%	3,9%
961,3 ug/dL	2,9%	4,3%

b) Detection limit: the detection limit is 15 ug/dL (9 umol/L) of ammonia.

c) Linearity: the reaction is linear up to 1700 ug/dL (1000 umol/L) of ammonia.

PARAMETERS FOR AUTOMATIC ANALYZERS

For programming instructions refer to the user manual of the analyzer in use.

WIENER LAB PROVIDES

- 1 x 20 ml (Cat. N° 1009693)

- 1 x 20 ml (Cat. N° 1009911)*

REFERENCES

- Tietz, N.W., Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, PA: W.B. Saunders Co; 1995:44.
- Young, D.S., Effects of preanalytical variables on clinical laboratory test. AACC Press. Third ed., 2007.
- Mondzac A, Ehrlich GE, Seegmiller JE - An enzymatic determination of ammonia in biological fluids - J. Lab. Clin. Med. 66/3:526-31, 1965.

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Wiener Laboratorios S.A.I.C.