Apo A-I

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(Cód. 1009947)*

REFERÊNCIAS

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.

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Apo A-I

Immunoturbidimetric method for quantitative determination of apolipoprotein A-I in serum

SUMMARY

Apolipoprotein A-I (Apo A-I) is the main component of high density lipoprotein (HDL). It is synthesized in the small intestine, liver and in small proportion in the kidneys and gonads. Apo A-I has an important physiological role: it acts as cofactor in lecithin cholesterol acyltransferase (LCAT) and also removes free cholesterol from cells.

The process of esterification of cholesterol catalyzed by LCAT is important in the reverse transport of cholesterol to the liver in order to be metabolized and excreted.

Additionally, Apo A-I has a role of prostacyclin stabilization and thus enhances vasodilation and inhibits platelet aggregation. These functions are important in protection against atherosclerosis.

Apo A-I levels are increased in the familial hyper-alpha lipoproteinemia during estrogen therapy, pregnancy and liver disease. Increased levels of Apo A-I are associated with a decreased risk of atherosclerosis.

Apo A-I levels are decreased in familial hypo-alpha-lipoproteinemia, in Tangier disease, cholestasis and sepsis.

PRINCIPLE

The apolipoprotein A-I reacts to the specific antibody forming insoluble immune complexes. The turbidity caused by these immune complexes is proportional to the apolipoprotein A-I concentration in the sample and may be spectrophotometrically measured.

PROVIDED REAGENTS

A. Reagent A: phosphate buffer, pH 7.4.

B. Reagent B: polyclonal antibodies anti-human apolipoprotein A-I (goat) in phosphate buffer, pH 7.4.

NON-PROVIDED REAGENTS

- Saline solution
- Wiener lab.'s **Apo Calibrator Turbitest AA**

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

The reagents are for "in vitro" diagnostic use.

All patient samples should be handled as capable of transmitting infectious diseases.

Use the reagents according to the working procedures for clinical laboratories.

Reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date stated on the box. Do not freeze.

SAMPLE

Serum

a) Collection: obtain in the usual way.

b) Additives: not required.

c) Known interfering substances: do not use hemolyzed, lipemic or contaminated samples. Before testing, particles in samples should be removed by centrifugation. No interferences are observed with hemoglobin up to 1000 mg/dl, triglycerides up to 2500 mg/dl, Bilirubin up to 20 mg/dl, heparin up to 50 mg/dl and sodium citrate up to 1000 mg/dl.

See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: sample should be preferably fresh. In case it cannot be processed immediately, the sample can be kept for up to 48 hours at 2-10°C or for longer period store at -20°C.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer.
- Square spectrophotometric cuvettes.
- Micropipettes and pipettes for measuring the stated volumes.
- Kahn or hemolysis tubes.
- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 340 nm
- Reaction temperature: room temperature (25°C). Temperature control is not critical, it can range between 22 and 30°C.
- Reaction time: 15 minutes
- Sample volume: 30 ul
- Final reaction volume: 1830 ul

Sample and reagents volumes may be proportionally changed without affecting the calculation factors.

PROCEDURE

CALIBRATION CURVE

In Kahn tubes dilute the Apo Calibrator Turbitest AA with saline solution 1:10, 1:20, 1:40, 1:80 and 1:160, using saline solution as the zero point.

Diluted Apo Calibrator	30 ul
Reagent A	1500 ul

Homogenize and measure the absorbance of each dilution at 340 nm (OD₁), setting the instrument to zero with distilled water. Then, add:

Reagent B	300 ul
------------------	--------

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD₂), setting the instrument to zero with distilled water.

Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each Apo Calibrator dilution, including the zero point. Draw on graph paper the ΔA absorbance differences based on the Apo Calibrator concentration in mg/dl (g/l).

SAMPLES PROCEDURE

Dilute the Samples 1:10 with saline solution.

Diluted Sample	30 ul
-----------------------	-------

Reagent A	1500 ul
------------------	---------

Homogenize and measure the absorbance at 340 nm (OD₁), setting the instrument to zero with distilled water. Then add:

Reagent B	300 ul
------------------	--------

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD₂), setting the instrument to zero with distilled water.

Intra-assay precision

Level	S.D.	C.V.
131.5 mg/dl	± 1.7 mg/dl	1.3%
281.2 mg/dl	± 8.9 mg/dl	3.2%

b) Detection limit: 20 mg/dl

c) Measuring range: 20 - 300 mg/dl

d) Prozone effect: not noted until 900 mg/dl Apolipoprotein A-I.

WIENER LAB. PROVIDES

60 ml: - 1 x 50 ml Reagent A
 - 1 x 10 ml Reagent B
 (Cat. N° 1009361)

60 ml: - 1 x 50 ml Reagent A
 - 1 x 10 ml Reagent B
 (Cat. N° 1009639)

60 ml: - 1 x 50 ml Reagent A
 - 1 x 10 ml Reagent B
 (Cat. N° 1009947)*

REFERENCES

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - «Effects of Drugs on Clinical Laboratory Tests», AACC Press, 5th ed., 2000.

CALCULATIONS

Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each sample tested. Interpolate this ΔA in the calibration curve to determine the concentration in mg/dl (g/l) corresponding to the sample under study. Samples with an absorbance above that of the Apo Calibrator Turbitest AA must be diluted with saline solution and processed again. Multiply the obtained result by the dilution factor.

QUALITY CONTROL METHOD

It is recommended the use of **Lipid Control**, separately provided by Wiener lab. The control should be processed in the same manner as samples.

REFERENCE VALUES

Men: 111 - 157 mg/dl (1.11 - 1.57 g/l)

Women: 126 - 182 mg/dl (1.26 - 1.82 g/l)

Each laboratory should set its own reference values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

It is recommended to perform a complete recalibration when changing reagent lot or when suggested by Quality Control. Avoid contamination to preserve the integrity of the reagents. Only use thoroughly clean and dry micropipettes for measurement.

PERFORMANCE

a) Reproducibility: replicates of samples containing different apolipoprotein A-I levels were assayed and the following results were obtained:

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Wiener Laboratorios S.A.I.C.