# ASO

## Calibrador

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# ASO

## Calibrador

For the calibration of anti-streptolysin O (ASO) determination by immunoturbidimetric method

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Nr kat. 1913263

### USES

**ASO Calibrador** has been designed to be used with Wiener lab's immunoturbidimetric reagents. Since it is lot-specific, check the given value on the label. That value has been standardized with the World Health Organization (WHO) standard.

### PROVIDED REAGENTS

**Calibrator:** serum with high concentration of anti-streptolysin O diluted in buffered saline solution and appropriate preservatives.

### INSTRUCTIONS FOR USE

**Calibrator:** ready to use.

### WARNINGS

The reagent is for "in vitro" diagnostic use.

Avoid ingestion and contact with the skin.

This calibrator has been prepared from non-reactive material to HBsAg, HCV and HIV. However, it should be handled as if it were an infective sample.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

### STABILITY AND STORAGE INSTRUCTIONS

**Calibrator:** stable at 2-10°C until the expiration date stated on the box.

Once opened is stable for up to 30 days at 2-10°C.

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Wiener Laboratorios S.A.I.C.