



# Bilirrubina Directa

AA

DPD method for direct bilirubin determination in serum or plasma

## SUMMARY

Bilirubin, a compound produced by the degradation of hemoglobin, is captured by the liver for its conjugation and bile excretion. Hepatocellular disorders or biliary obstructions may lead to hyperbilirubinemia. Erythroblastosis fetalis or hemolytic anemia of the newborn is a pathology produced by mother-fetus blood incompatibility, in which an excessive destruction of red blood cells occurs. This results in a severe increase of serum bilirubin, with the consequent risk of pigment diffusion to the central nervous system. Therefore, bilirubin determination in these newborns is extremely important.

## PRINCIPLE

Direct bilirubin reacts with dichlorophenyldiazonium salt (DPD) forming a red azocompound in an acid solution.

## PROVIDED REAGENTS

- A. Reagent A:** vials containing dichlorophenyldiazonium salt.  
**B. Reagent B:** aqueous solution containing 90 mmol/l sulfamic acid.  
**C. Reagent C** (Reagent for Sample Blank): aqueous solution containing 90 mmol/l sulfamic acid.

## Final concentrations

sulfamic acid ..... 90 mmol/l  
sodium chloride ..... 6.6 g/l  
dichlorophenyldiazonium ..... 0.1 mmol/l

## NON-PROVIDED REAGENTS

Wiener lab.'s **Calibrador A plus**.

## INSTRUCTIONS FOR USE

**Reagent B and C:** ready to use.

**Working Reagent:** reconstitute each vial of Reagent A with 10 ml of Reagent B. Cap and shake until complete dissolution. Homogenize and date.

## WARNINGS

Reagents are for "in vitro" diagnostic use.  
Reagent B and Reagent C are irritant. H315+H320: Causes skin and eye irritation. P262: Do not get in eyes, on skin, or on clothing. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P280: Wear protective gloves/protective clothing/eye protection/face protection. Use the reagents according to the working procedures for clinical laboratories. The reagents and samples should be discarded according to the local regulations in force.

## INSTABILITY OR DETERIORATION OF REAGENTS

Absorbance readings of the Blank Reagent, higher than 0.100 O.D. (at 546 nm) indicate its deterioration. Discard in such case.

Turbidity indicates reagents deterioration. Discard in such case.

## STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable in refrigerator (2-10°C) until the expiration date shown on the box.

**Working Reagent:** stable for 21 days in refrigerator (2-10°C).

## SAMPLE

Serum or plasma

**a) Collection:** obtain as usual. Protect from natural or artificial light covering the tube with black paper.

**b) Additives:** if plasma is used, use heparin as collection.

**c) Known interfering substances:** no interference have been observed by hemoglobin up to 180 mg/dl, triglycerides up to 550 mg/dl (5,50 g/l) or heparin up to 20 U/ml.

See Young, D.S. in References for effect of drugs on the present method.

**d) Stability and storage instructions:** sample should be preferably fresh. If assay is not performed immediately, serum can be stored up to 48 hours in refrigerator (2-10°C) and whole blood no more than 24 hours in refrigerator (2-10°C) or 12 hours at room temperature (< 25°C).

The action of light is capable of destroying up to a 50% of the bilirubin present in the sample. Consequently, it should be carefully protected from light.

## REQUIRED MATERIAL (non-provided)

- Spectrophotometer or photocolormeter
- Micropipettes and pipettes for measuring the stated volumes
- Tubes or spectrophotometric cuvettes
- Stopwatch

## ASSAY CONDITIONS

- Wavelength: 546 nm in spectrophotometer or 520-550 nm in photocolormeter with green filter.
- Reaction temperature: room temperature (< 25°C)
- Reaction time: 10 minutes
- Sample volume: 40 ul
- Final reaction volume: 0.54 ml

## PROCEDURE

In 4 photocolormeter tubes labeled B<sub>C</sub> (Calibrator Blank),

C (Calibrator), B<sub>U</sub> (Unknown Blank) and U (Unknown), place:

	B	C	B <sub>U</sub>	U
<b>Reagent C</b>	0.5 ml	-	0.5 ml	-
<b>Reagent B</b>	-	0.4 ml	-	0.4 ml
<b>Calibrator</b>	40 ul	40 ul	-	-
<b>Sample</b>	-	-	40 ul	40 ul
<b>Working Reagent</b>	-	0.1 ml	-	0.1 ml

Mix and incubate for 10 minutes at room temperature (< 25°C). Read in spectrophotometer at 546 nm or in photocolorimeter with green filter (520-550 nm), setting the instrument to zero with the Reagent Blank (B<sub>R</sub>). Volumes may be proportionally increased.

Note: process a Reagent Blank (B<sub>R</sub>) for each series of determinations mixing 0.4 ml Reagent B + 40 ul distilled water + 0.1 ml Working Reagent.

### STABILITY OF FINAL REACTION

Absorbance should be read at exactly 10 minutes.

### CALCULATIONS

Direct Bilirubin (mg/l) = (U - B<sub>U</sub>) x f

$$f = \frac{X^* \text{ mg/l}}{C - B_c}$$

\*direct bilirubin concentration in Wiener lab.'s **Calibrador A plus**

Example:

B<sub>C</sub> = 0.008

C = 0.135

B<sub>U</sub> = 0.009

U = 0.029

X = 23.9 mg/l

$$f = \frac{23.9 \text{ mg/l}}{0.135 - 0.008} = 188.2 \text{ mg/l}$$

The calculated colorimetric factor (f) may also be employed, using Wiener lab's **Bilirrubina Standard** with Wiener lab.'s **Bilirrubina Total AA** (DPD method) technique.

### QUALITY CONTROL METHOD

Each time the test is performed, analyze two levels of a quality control material (**Standatrol S-E 2 niveles**) with known direct bilirubin concentration.

### REFERENCE VALUES

Direct bilirubin in serum or plasma:

Adults: up to 2 mg/l

It is recommended that each laboratory establishes its own reference values.

### UNITS CONVERSION

Bilirubin (mg/l) = Bilirubin (mg/dl) x 10

Bilirubin (mg/dl) x 17.1 = Bilirubin (umol/l)

### PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

The action of light, both on sample and on the standard solutions, are capable of destroying up to 50% of the bilirubin in one hour.

### PERFORMANCE

The assays were performed in an Express Plus analyzer<sup>(\*)</sup>. If using by manual procedure, user must validate that similar performance to that stated below is obtained.

**a) Reproducibility:** CLSI's (ex NCCLS) protocol EP15-A was applied. The following results were obtained:

#### Intra-assay precision

Level	S.D.	C.V.
1.4mg/l	± 0.05 mg/l	4.03 %
18.0 mg/l	± 0.24 mg/l	1.33 %

#### Total precision

Level	S.D.	C.V.
1.4 mg/l	± 0.06 mg/l	5.09 %
18.0 mg/l	± 0.36 mg/l	2.00 %

**b) Linearity:** applying the EP6-P protocol, a linearity of 140 mg/l was obtained. For higher values, dilute the sample 1:2 or 1:4 with saline solution, repeat the determination and multiply the obtained result by 2 or 4 accordingly.

**c) Analytical sensitivity:** 1.1 mg/l.

### PARAMETERS FOR AUTOANALYZERS

For programming instructions check the user's manual of the autoanalyzer in use.

For calibration, it must be used Wiener lab.'s **Calibrador A plus**.

The kit provides the Reagent for the Sample Blank (Reagent C) which is only required in some autoanalyzers. In other analyzers this Reagent is not necessary. For further information consult your analyzer's application manual for Wiener lab. reagents.

### WIENER LAB. PROVIDES

- 200 ml (Cat. 1009306): - 4 vials Reagent A  
 - 4 x 10 ml Reagent B  
 - 4 x 40 ml Reagent C

- 200 ml (Cat. 1120006): - 4 vials Reagent A  
 - 4 x 50 ml Reagent B  
 - 2 x 100 ml Reagent C

### REFERENCES

- Bartels, H. and Boehmer, M. - Z. Klin. Chem. Klin. Biochem. (1971).
- Colombo, J.; Peheim, E.; Kyburz, S. and Hoffman - Clin. Chim. Acta 51/2:217 (1974).
- Botwell, J. H. - Clin. Chem. 10/3:197 (1964).
- Watson, D. - Clin. Chem. 7/6:603 (1961).

- Ichida, T. and Nobuoka, M. - Clin. Chim. Acta 19/2:249 (1968).
- Rand and Di Pascua - Clin. Chem. 8/6 (1962).
- Zaroda, R. - Am. J. Clin. Path. 45/1:70 (1996).
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 3rd ed. (1990).
- Tietz Textbook of Clinical Chemistry - Saunders Co., 3<sup>rd</sup> ed. (1999).
- CLSI: Clinical and Laboratory Standards Institute (ex-NCCLS). Document "Evaluation of the Linearity of Quantitative Analytical Methods", EP6-P (1986).
- CLSI: Clinical and Laboratory Standards Institute (ex-NCCLS). Document "Evaluation of Precision Performance of Clinical Laboratory Devices", EP5-A (1999).

## SYMBOLS

The following symbols are used in the packaging for Wiener lab. diagnostic reagents kits.

 This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

 Authorized representative in the European Community

 "In vitro" diagnostic medical device

 Contains sufficient for <n> tests

 Use by

 Temperature limitation (store at)

 Do not freeze

 Biological risks

 Volume after reconstitution

 Contents

 Batch code

 Manufactured by:

 Harmful

 Corrosive / Caustic

 Irritant

 Consult instructions for use

 Calibrator

 Control

 Positive Control

 Negative Control

 Catalog number

 Wiener Laboratorios S.A.I.C.  
Riobamba 2944  
2000 - Rosario - Argentina  
<http://www.wiener-lab.com.ar>  
Dir. Téc.: Viviana E. Cétola  
Biochemist  
A.N.M.A.T. Registered product  
Cert. N°: 4100/00

 **Wiener lab.**

2000 Rosario - Argentina