



Calibrador A *plus*

Clinical chemistry calibrator for autoanalyzers

Kit lot: 1508169560
Calibrador A plus lot: 169560
Reagent A lot: 169560

USES

The calibrator has been designed to be used with Wiener lab.'s reagents in autoanalyzers.

Check the value chart for the constituents, since they are lot specific.

PRINCIPLE

Calibrador A plus contains the compounds usually determined in clinical chemistry laboratories.

It must be noted that the assigned values for the calibrator's different constituents have been obtained with Wiener lab.'s methods and reagents, thus, the obtained results will only be comparable to the list if the corresponding methods and reagents have been used.

PROVIDED REAGENTS

Calibrator: lyophilized serum containing metabolites in the appropriate concentrations to ensure a correct calibration of autoanalyzers.

A. Reagent A: 25 mmol/l sodium carbonate solution, pH 10.

INSTRUCTIONS FOR USE

- Open the vial, slowly pulling out the rubber stopper to avoid any loss of the lyophilized material.
- Add exactly 3 ml Reagent A (using a burette or a pipette with double calibration mark).
- Cap and mix by gentle inversion, avoiding the formation of foam. Do not shake.
- Let dissolve for 30 minutes at room temperature, mixing by inversion from time to time.
- Mix by inversion immediately before use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

This calibrator has been prepared from non-reactive material to HBsAg, HCV and HIV. However, it should be handled as if it were an infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

The **Provided Reagents** are stable in refrigerator (2-10°C) until the expiration date shown on the box.

Reconstituted Calibrator: is stable for 8 hours at room temperature (<25°C), 2 days refrigerated (2-10°C) or 30 days frozen (-20°C).

In the dark, bilirubin is stable 4 hours at room temperature (<25°C), 8 hours refrigerated (2-10°C) or 2 weeks frozen (-20°C).

Avoid repeated freezing and thawing.

INSTABILITY OR DETERIORATION OF REAGENTS

Any variation of the organoleptic features of the Calibrator may be a sign of its deterioration.

PROCEDURE

Reconstituted Calibrator should be used in the same manner as an unknown sample, following the autoanalyzer parameters that come with the reagent kit in use.

REFERENCE VALUES

Calibrator values have been established using Reference Material from the National Bureau of Standards, when it was possible.

PROCEDURE LIMITATIONS

Mistakes in the reconstitution or conservation may cause erroneous results.

See PROCEDURE LIMITATIONS on the package insert of the kit in use.

WIENER LAB. PROVIDES

- 2 x 3 ml (Cat. 1918005).

CALIBRADOR A PLUS METABOLITES ASSIGNED VALUES (37°C)

Analyte	Method	Value	Unit
Uric acid	Uricostat enzimático AA	4.74	mg/dL
	Uricostat enzimático AA líquida	4.74	mg/dL
Albumin	Albúmina AA	3.55	g/dL
Direct bilirubin	Bilirrubina	2.57	mg/dL
	Bilirrubina Directa AA	2.42	mg/dL
	Bilirrubina Directa AA líquida	2.19	mg/dL
Total bilirubin	Bilirrubina	5.11	mg/dL
	Bilirrubina Total AA	4.28	mg/dL
	Bilirrubina Total AA líquida	4.27	mg/dL
Calcium	Ca-Color AA	10.6	mg/dL
	Ca-Color Arsenazo III AA	10.6	mg/dL
Cholesterol	Colestat enzimático AA	160	mg/dL
	Colestat enzimático AA líquida	159	mg/dL
Creatinine	Creatinina cinética AA líquida - Conventional technique	3.69	mg/dL
	Creatinina cinética AA líquida - Compensate technique	3.99	mg/dL
	Creatinina enzimática AA líquida	3.80	mg/dL
Phosphorus	Fosfatemia UV AA	5.81	mg/dL
Glucose	Glicemia enzimática AA	187	mg/dL
	Glicemia enzimática AA líquida	188	mg/dL
HDL Cholesterol	HDL Colesterol Reactivo Precipitante	176	mg/dL
	HDL Colesterol FT	560	mg/dL
Iron	Fer-color AA	201	ug/dL
	Fer-color AA líquida	189	ug/dL
Lactate	Lactate	29.3	mg/dL
LDL Cholesterol	LDL Colesterol Reactivo Precipitante	240	mg/dL
Magnesium	Mg-Color AA	2.73	mg/dL
	Magnesium CPZ	2.73	mg/dL
Total proteins	Proteínas Totales AA	5.35	g/dL
	Total Protein	5.11	g/dL
Triglycerides	TG Color GPO/PAP AA	131	mg/dL
	TG Color GPO/PAP AA líquida	130	mg/dL
UIBC	UIBC/TIBC AA líquida	146	ug/dL
Urea	Urea UV cinética AA	101	mg/dL
	Urea UV cinética AA líquida	101	mg/dL

CALIBRADOR A PLUS ENZYMES ASSIGNED VALUES (37°C)

Analyte	Method	Value	Unit
Alanine aminotransferase (GPT/ALT)	GPT (ALT) UV AA	106	U/L
	GPT (ALT) UV AA líquida	106	U/L
Amylase	Amilasa 405 cinética AA	263	U/L
	Amilasa 405 cinética AA líquida	258	U/L
Aspartato aminotransferase (GOT/AST)	GOT (AST) UV AA	109	U/L
	GOT (AST) UV AA líquida	109	U/L
Cholinesterase	Colinesterasa AA	5134	U/L
	Cholinesterase	4804	U/L
Creatine kinase	CK-NAC UV AA	338	U/L
	CK-NAC UV AA líquida	351	U/L
Creatine kinase - MB	CK-MB NAC UV AA líquida	312	U/L
Acid phosphatase	Fosfatasa Acida Total y Prostática cinética - FANP	9.71	U/L
	Fosfatasa Acida Total y Prostática cinética - FAT	25.7	U/L
Alkaline phosphatase	ALP 405 AA líquida	466	U/L
γ -Glutamyl transferase	γ -G test cinética AA	100	U/L
	γ -G test cinética AA líquida	101	U/L
Lactate deshydrogenase	LDH-P UV AA	523	U/L
	LDH-P UV AA líquida	515	U/L
	LDH-L	257	U/L
Lipase	Lipasa AA líquida	99	U/L

Symbols


The following symbols are used in the packaging of Wiener lab. diagnostic reagent kits.



This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices



Manufactured by:



Authorized representative in the European Community



Harmful



"In vitro" diagnostic medical device



Corrosive / Caustic



Contains sufficient for <n> tests



Irritant



Use by



Consult instructions for use



Temperature limitation (store at)



Do not freeze



Calibrator



Biological risks



Control



Volume after reconstitution



Positive Control



Contents




Negative Control



Batch code



Catalog number

 Wiener Laboratorios S.A.I.C.
 Riobamba 2944
 2000 - Rosario - Argentina
<http://www.wiener-lab.com.ar>
 Tech.Dir.: Viviana E. Cétola
 Biochemist
 A.N.M.A.T. Registered product
 Cert. N°: 1976/97



Wiener lab.

2000 Rosario - Argentina