# Calibrador Proteínas

*nivel alto*

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Nr kat. 1913261  
Nr kat. 1913268

### USES

**Calibrador Proteínas *nivel alto*** has been designed to be used with Wiener lab's turbidimetric reagents. Since they are lot specific, check the value chart given for the constituents.

### PROVIDED REAGENTS

**Calibrador Proteínas *nivel alto***: inactivated human serum with the appropriate preservatives.

### INSTRUCTIONS FOR USE

**Calibrador Proteínas *nivel alto***: ready to use.

### WARNINGS

The reagent is for "in vitro" diagnostic use.

Avoid ingestion and skin contact.

This calibrator has been prepared from non-reactive material to HBsAg, HCV and HIV. However, it should be handled as if it were an infective sample.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

### STABILITY AND STORAGE INSTRUCTIONS

**Calibrador Proteínas *nivel alto***: stable at 2-10°C until the expiration date stated on the box. Once opened is stable for up to 30 days at 2-10°C. Do not freeze.

### ASSIGNED VALUES

The values are based on the International Federation of Clinical Chemistry (IFCC) RPPHS/CRM 470 reference preparation.

### WIENER LAB. PROVIDES

- 1 x 1 ml (Cat. 1913261).
- 2 x 1 ml (Cat. 1913268).

### REFERENCES

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Tietz Textbook of Clinical Chemistry - Burtis, C.; Ashwood, E. 5th Edition WB Saunders, 2001.
- Hafner G et al. - Effects of standardization with the new international reference preparation for proteins in human serum on method comparability and reference values. - Clin. Lab. 41:743-748, 1995.
- EP5-A Vol. 24 N° 25 del CLSI - Evaluation of precision performance of quantitative measurement methods (approved guideline - second edition).

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Wiener Laboratorios S.A.I.C.