# CK-MB

## Control - 3 niveles

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# CK-MB

## Control - 3 niveles

For precision control in the MB isoenzyme determination  
of Creatine Kinase

### USES

**CK-MB Control 3 niveles** is a lyophilized control based on bovine serum albumin using human CK-MB isoenzyme and is designed to be used in the precision control of **CK-MB NAC UV AA líquida** and **CK-MB DS UV unitest** kits.

### PROVIDED REAGENTS

**CK-MB Control level 1:** lyophilized control based on bovine serum albumin with the addition of human CK-MB isoenzyme in low and/or close to normal activity levels.

**CK-MB Control level 2:** lyophilized control based on bovine serum albumin with the addition of human CK-MB isoenzyme close to normal activity levels.

**CK-MB Control level 3:** lyophilized control based on bovine serum albumin with the addition of human CK-MB isoenzyme in pathological activity levels.

### NON-PROVIDED REAGENTS

Distilled or demineralized water.

### INSTRUCTIONS FOR USE

Carefully open the vial removing the rubber cap to avoid the loss of lyophilized material. Reconstitute adding the quantity of distilled or demineralized water stated on the label. Cap and homogenize gently by inversion until complete dissolution of the material. Avoid foam formation.

Add the necessary volume to a sample container and test as a patient sample. It is advisable to process the control daily, in parallel to patient samples, and to establish control intervals according to the needs of each particular laboratory. The results should be within the defined intervals. Each laboratory should establish the guidelines for corrective measures if the values are established outside the interval. Each laboratory should establish quality control procedures to meet current local regulations or their accreditation requirements.

### WARNINGS

Reagents are for "in vitro" diagnostic use.

The Controls have been prepared from non-reactive material to Hepatitis B surface antigen (HBsAg), antibodies to Hepatitis C virus (HCV) and Human Immunodeficiency Virus (HIV). However, because no known test method can offer complete assurance that infectious agents are absent, the controls and the samples should be handled as potentially infectious material.

Use gloves and adequate eye protection.

Avoid ingestion and/or direct contact with the skin and mu-

cos membranes. In such case, proceed according to the instructions of the competent health authorities.

Use the reagents following the usual work precautions at the clinical laboratory.

All reagents and samples should be discarded according to the current local regulations.

### STABILITY AND STORAGE INSTRUCTIONS

The lyophilized controls are stable at 2-10°C until the expiration date stated on the box. Once reconstituted, they are stable for up to 3 days at 2-10°C or for up to 2 days at 25°C. After use, closed vials firmly and keep them at 2-10°C and protected from light. When the reconstituted controls are stored at -20°C, they are stable for 3 months in such conditions. Aliquot and freeze them immediately after reconstitution. Avoid repeated freezing and thawing.

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Wiener Laboratorios S.A.I.C.