



# CK-NAC UV AA

unitest y AA


Optimized UV method (IFCC) for the determination of Creatine Kinase (CK) in serum or plasma

## SUMMARY

Creatine Kinase is an intracellular enzyme mostly found in skeletal and heart muscles, as well as in the brain. Therefore, an increase in serum activity indicates cellular damage. In Acute Myocardial Infarction (AMI), the increase of CK serum activity starts 2 to 6 hours after the onset and reaches its peak after 18 to 24 hours. Since peaks may be 20 times the highest peak increase, this is perhaps the most sensitive test for AMI diagnosis.

## PRINCIPLE

The reaction system is as follows:



In the reaction system, NAC works as activator of the Creatine Kinase, recommended by the IFCC.

## PROVIDED REAGENTS

**Reagent B:** Imidazole buffer solution, pH 6.7.

**Reagent A:** vial containing dry substrate in the required quantities for the following final concentrations:

Imidazole .....	100 mmol/l; pH 6.7
Creatine phosphate .....	30 mmol/l
ADP .....	2 mmol/l
Glucose .....	20 mmol/l
NADP .....	2 mmol/l
Hexokinase .....	≥ 2500 U/l
Glucose-6-phosphate dehydrogenase .....	≥ 2000 U/l
Magnesium acetate .....	10 mmol/l
AMP .....	5 mmol/l
Di (adenosine-5') pentaphosphate .....	10 umol/l
N-Acetyl-cysteine .....	20 mmol/l

## INSTRUCTIONS FOR USE

**Reagent B:** ready to use.

**Reagent A:** reconstitute according to the following:

- CK-NAC UV AA: dissolve each Reagent A vial with 20 ml Reagent B. Date.
- CK-NAC UV unitest: dissolve each Reagent A vial with 2.5 ml Reagent B. Date.

## WARNINGS

Reagents are for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable in refrigerator (2-10°C) until the expiration date shown on the box.

**Reconstituted Substrate:** stable in refrigerator for 20 days (2-10°C) or 3 days at room temperature from reconstitution date.

## INSTABILITY OR DETERIORATION OF REAGENTS

When the spectrophotometer has been set to zero with distilled water, absorbance readings of the reconstituted Reagent A higher than 0.800 O.D. (at 340 nm) indicate its deterioration.

## SAMPLE

Serum or plasma

**a) Collection:** obtain in the usual way.

**b) Additives:** when using plasma, use only heparin (concentration < 15 U/l) or EDTA as anticoagulant.

**c) Known interfering substances:** no interferences are observed from bilirubin up to 90 mg/l, triglycerides up to 13 g/l, heparin up to 50 IU/l. Sera with visible hemolysis (hemoglobin concentration > 008 g/dl) produce interferences.

See Young, D.S. in References for effect of drugs on the present method.

**d) Stability and storage instructions:** sample should be preferably fresh. It can be kept up to 1 week in refrigerator (2-10°C) without preservatives.

## REQUIRED MATERIAL (non-provided)

- Spectrophotometer.
- Micropipettes and pipettes for measuring the stated volumes
- Water bath at the temperature indicated under PROCEDURE.
- Stopwatch.

## ASSAY CONDITIONS

(Increase of Absorbance)

- Wavelength: 340 nm (Hg 334 or 366)
- Reaction temperature: 25, 30 or 37°C. See the REFERENCE VALUES corresponding to each temperature.
- Reaction time: varies according to the selected procedure.

## PROCEDURE

Set instruments to zero with distilled water.

### A) 30 - 37°C

I- MACROTECHNIQUE

In a cuvette at 30-37°C, place:

**Reconstituted Reagent A** 2.5 ml

Pre-incubate 3-4 minutes. Then add:

**Sample** 50 ul

Mix immediately and wait 3 minutes. Adjust absorbance to a reference reading and simultaneously start stopwatch. Record absorbance 1, 2 and 3 minutes after first reading. Determine average change in Absorbance/min ( $\Delta A/\text{min}$ ), subtracting each reading from the previous one and averaging these values. Use this mean for calculations.

**II- MICROTECHNIQUE**

In a cuvette at 30-37°C, place:

**Reconstituted Reagent A** 1 ml

Pre-incubate 3-4 minutes. Then add:

**Sample** 20 ul

Mix immediately and wait 3 minutes. Follow A-I procedure.

**B) 25°C**

**I- MACROTECHNIQUE**

Since activity at this temperature is lower, use 100 ul sample to reach required sensitivity. Start recording readings 4 minutes after sample addition. Follow A-I procedure.

**II- MICROTECHNIQUE**

Follow the steps described in A-II, but using 50 ul sample.

**CALCULATIONS**

$CK\ U/I = \Delta A/\text{min} \times \text{factor}$

In each case the corresponding calculation factor should be used, as shown on the table below:

Temperature \ Wavelength	30-37°C	25°C	
	I or II	I	II
340 nm	8,095	4,127	3,333
334 nm	8,252	4,207	3,398
366 nm	15,000	7,647	6,176

**QUALITY CONTROL METHOD**

Each time the test is performed, analyze two levels of a quality control material (**Standatrol S-E 2 niveles**) with known creatine kinase activity.

**REFERENCE VALUES**

Temperature	25°C	30°C	37°C*
Men	up to 80 U/I	up to 130 U/I	up to 195 U/I
Women	up to 70 U/I	up to 110 U/I	up to 170 U/I

\*Calculated

It is recommended that each laboratory establishes its own reference values.

**SI SYSTEM UNITS CONVERSION**

Creatine kinase (U/I) x 0.017 = Creatine kinase (ukat/l)

**PROCEDURE LIMITATIONS**

See Known Interfering Substances under SAMPLE.

**PERFORMANCE**

**a) Reproducibility:** simultaneously processing replicates of the same sample in the same day, were obtained:

Level	S.D.	C.V.
144.5 U/I	± 3.36 U/I	2.33 %
450.4 U/I	± 5.98 U/I	1.33 %

Performing the same assay on different days, were obtained:

Level	S.D.	C.V.
147.5 U/I	± 3.13 U/I	2.12 %
451.2 U/I	± 6.89 U/I	1.53 %

**b) Detection limit:** depends on the photometer used and the wavelength. In spectrophotometer with 1 cm optical length square cuvettes, for a  $\Delta A/\text{min}$  of 0.001, the minimum detectable activity change will be of 8 U/I (at 340 nm and 30 or 37°C).

**c) Dynamic range:** the useful reading range is extended up to 0.250 O.D. (340 nm) or 0.140 O.D. (366 nm). For higher values dilute the sample with saline solution, correcting the results accordingly.

**PARAMETERS FOR AUTOANALYZERS**

For programming instructions check the user manual of the autoanalyzer in use.

**WIENER LAB. PROVIDES**

**CK-NAC UV AA:**

- 3 x 20 ml (60 ml Buffer) (Cat. N° 1271303).
- 3 x 20 ml (3 x 20 ml Buffer) (Cat. N° 1009309).
- 10 x 20 ml (200 ml Buffer) (Cat. N° 1271353).

**CK-NAC UV unitest:**

- 20 vials x 2.5 ml Reagent (Cat. N° 1271351).

**REFERENCES**

- D.G.K.C. - Z. Klin. Chem. 10:281 (1972).
- S.S.C.C. - Scand. J. Clin. Lab. Invest. 33:291 (1974).
- I.F.C.C. - Clinica Chimica Acta 105:147 F (1980).
- I.F.C.C. - Ann. Biol. Clin. 44/4:419 (1986).
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

## SYMBOLS

The following symbols are used in packaging for Wiener lab. diagnostic reagents kits.



This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

	Authorized representative in the European Community
	"In vitro" diagnostic medical device
	Contains sufficient for <n> tests
	Use by
	Temperature limitation (store at)
	Do not freeze
	Biological risks
	Volume after reconstitution
	Contents
	Batch code
	Manufactured by:
	Harmful
	Corrosive / Caustic
	Irritant
	Consult instructions for use
	Calibrator
	Control
	Positive Control
	Negative Control
	Catalog number

 Wiener Laboratorios S.A.I.C.  
Riobamba 2944  
2000 - Rosario - Argentina  
<http://www.wiener-lab.com.ar>  
Dir. Téc.: Viviana E. Cétola  
Bioquímica  
A.N.M.A.T. Registered product  
CK-NAC UV unitest:  
Disp. N° 7448/87 - 370/00  
CK-NAC UV AA:  
Disp. N° 5999/83-5234/98-955/02



**Wiener lab.**

2000 Rosario - Argentina