CMD 80 Detergent

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- 1 x 2 L (Cód. 1009676)

CMD 80 Detergent

For cleaning probes, mixers and cuvettes in automated clinical chemistry analyzer

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Nr kat. 1009666
Nr kat. 1009676

USE

CMD 80 Detergent is intended to be used in automated clinical chemistry analyzers (Wiener lab. CMD 800 series) for removing proteins, lipids, ions and other chemistry reaction residues on probes, mixers and cuvette surfaces.

PROVIDED REAGENTS

CMD 80 Detergent: 5% Potassium hydroxide, 1% Citrate buffer and surfactants.

INSTRUCTIONS FOR USE

Ready to use.

WARNINGS

- Reagent is for "in vitro" diagnostic use.
- CMD 80 Detergent is corrosive. H315+H320: Causes skin and eye irritation. H314: Causes severe skin burns and eye damage. P262: Do not get in eyes, on skin, or on clothing. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P280: Wear protective gloves/protective clothing/eye protection/face protection.
- Use the reagents according to the working procedures for clinical laboratories.
- The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

CMD 80 Detergent: stable at room temperature (2-30°C) until de expiration date stated on the box.

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Wiener Laboratorios S.A.I.C.