



# Factor V

## *Deficient Plasma*

For coagulometric determination of Factor V with one-stage method

### SUMMARY

Factor V (FV) is an essential protein in the homeostatic regulation of blood given its important role as cofactor in the mechanisms of coagulation and anticoagulation. It is synthesized primarily in the liver as a single chain glycoprotein suffering numerous post translational modifications before being secreted into the blood. Besides circulating free in plasma, FV is present in  $\alpha$  granules of platelets, thus representing 25% of the circulating FV. During coagulation, platelet FV is secreted as a result of platelet activation.

FV deficiencies can be inherited or acquired. The inherited form is autosomal recessive and its phenotypic expression is variable. Heterozygotes being generally asymptomatic and homozygotes showing symptoms of varying intensity, from mild to severe. The acquired deficit may be due to the presence of specific inhibitors or be associated with other factors such as liver damage, disseminated intravascular coagulation, etc.

One of the most common inherited causes of venous thrombosis is a mutation of FV known as FV Leiden, which is manifested by an Activated Protein C Resistance (APCR) phenotype.

### PRINCIPLE

Quantitative determination of FV consists of measuring the clotting time of a diluted sample containing the factor to be determined, using a deficient plasma providing the other factors in adequate levels, except FV, in the presence of calcium thromboplastin (one-step Prothrombin Time). The clotting time obtained is inversely proportional to the activity of FV present in the sample.

Factor V Deficient Plasma can be used in the predilution of patient samples in determining APCR.

This method can be used with any instrument capable of performing factor testing based on Prothrombin Time.

### PROVIDED REAGENTS

**A. Reagent A:** lyophilized human plasma deficient in factor V obtained by immunoabsorption with a coagulation activity <1% of FV.

### INSTRUCTIONS FOR USE

Dissolve Reagent A in the distilled water volume stated on the label. Let stand for 30 minutes at room temperature and then homogenize the solution by gentle shaking before use.

### NON-PROVIDED REAGENTS

- Distilled water

- Imidazole Buffer from Wiener lab.
- Soluplastin from Wiener lab or other PT reagents.
- Coagulation Control N and Coagulation Control P from Wiener lab.
- Coagulation Calibrator from Wiener lab.

### WARNINGS

Reagents are for "in vitro" diagnostic use.

**Factor V Deficient Plasma** has been prepared from non-reactive material for HBsAg, HCV and HIV. However, like blood samples, it should be handled as if they were infectious material.

Use the reagents keeping the usual work precautions in clinical laboratories.

All reagents and samples should be discarded according to local regulations.

### STABILITY AND STORAGE INSTRUCTIONS

**Factor V Deficient Plasma** is stable in refrigerator (2-10°C) until the expiration date shown on the box.

Once reconstituted the reagent is stable for 2 hours at room temperature (<25°C) or for 1 month frozen (-20°C). Avoid repeated freezing and thawing.

The frozen reagent should be thawed for at least 10 minutes at 37°C and homogenized before use.

### SAMPLE

Citrated plasma

**a) Collection:** obtain blood carefully (avoiding stasis or trauma) and place in a tube with anticoagulant in 9 + 1 exact proportion (example: 4.5 ml blood + 0.5 ml anticoagulant). Mix gently. Centrifuge for 15 minutes and separate the plasma within 30 minutes. It is advisable to collect with plastic syringes.

**b) Additives:** to obtain plasma Anticoagulant TP from Wiener lab should be used or 130 mmol/l (3.8%) or 109 mmol/l (3.2%) sodium citrate.

**c) Known Interfering Substances:**

- Avoid using EDTA or heparin for plasma collection.
- Contamination, visible or not, are the cause of falsely prolonged times.
- Hemolysis and visible lipemias obstruct photo-optical measurement of results.

Refer to the literature of Young for the effects of drugs on the present method.

**d) Stability and storage instructions:** plasma should be kept at room temperature until testing. This period should not last longer than 4 hours. In case plasma could not be

processed at once, it can be frozen up to 2 weeks at -20 °C. In this case the sample should be frozen immediately and should be thawed quickly at 37°C, no longer than 10 minutes.

**REQUIRED MATERIAL** (non-provided)

- Hemolysis tubes.
- Pipettes or micropipettes for measuring stated volumes.
- Stopwatch.
- Water bath at 37°C.
- Semiautomatic or automatic coagulometer.

**PROCEDURE**

**I- CALIBRATION CURVE**

- 1- Dilute 1/5 Coagulation Calibrator from Wiener lab in Imidazole Buffer (1 part Calibrator + 4 parts Imidazole Buffer) and prepare a series of geometric dilutions in Imidazole Buffer (1: 1 to 1:32). The 1: 1 ratio corresponds to 1/5 calibrator dilution, therefore dilutions shall be: 1/5; 1/10, 1/20, 1/40, 1/80 and 1/160. Calibrator diluted 1:10 represents 100% of the value assigned to FV in the table of calibrator values.
- 2- Preheat PT reagent (Soluplastin from Wiener lab) at 37°C.
- 3- In a hemolysis tube place:

<b>Reagent A</b>	100 ul
<b>Calibrator dilutions</b>	100 ul

- 4-Mix and incubate for 1 minute at 37°C
- 5- Start the stopwatch with the addition of 200 uL pre-heated PT Reagent and record the clotting time.
- 6- Calculate the average clotting time for each dilution, in duplicate.
- 7- Plot the calibration curve clotting times depending on FV activity on log-log paper. Join most of the points represented with a straight line. The final line should contain at least 3 consecutive points.

FV activity at each dilution curve is determined by multiplying the Coagulation Calibrator FV% by the dilution factor. Example for an assigned value of 98% to Coagulation Calibrator:

Co- agulation Calibrator Dilution (prediluted 1/5)	Final Dilution	FV Cal- culation (%)	Factor V (%)
1:1	1:5	98 x 2	196
1:2	1:10	98 x 1	98
1:4	1:20	98 x 0.5	49
1:8	1:40	98 x 0.25	24.5
1:16	1:80	98 x 0.125	12.3
1:32	1:160	98 x 0.063	6.1

**II- PATIENT SAMPLES**

- 1- Prepare 1:10 dilutions from patients plasmas in Imidazole Buffer (1 part sample + 9 parts Imidazole Buffer).
- 2- Preheat PT reagent (Soluplastin from Wiener lab) at 37°C.
- 3- In a hemolysis tube place:

<b>Factor V Deficient Plasma</b>	100 ul
<b>Diluted sample</b>	100 ul

- 4- Mix and incubate for 1 minute at 37°C
- 5- Start the stopwatch with the addition of 200 uL pre-heated PT Reagent and record the clotting time.
- 6- Repeat the determination and average the results for each sample.

**CALCULATIONS**

The values of the 1:10 diluted plasma samples are interpolated on the calibration curve.

**INTERPRETATION OF RESULTS**

If FV activity obtained by direct interpolation of the calibration curve is below the lowest point of the curve, repeat the determination of the sample with a lower dilution (1: 5), multiplying the result by 0.5.

If the activity of FV obtained by direct interpolation of the calibration curve is higher than the highest point of the curve, repeat the determination of the sample with a higher dilution (1:20), multiplying the result by 2.

**QUALITY CONTROL METHOD**

Coagulation Control N and Coagulation Control P from Wiener lab.  
The control shall be processed in the same manner as the samples.

**REFERENCE VALUES**

70-120%

Each laboratory should establish its own reference values from the techniques and instruments used.

**PROCEDURE LIMITATIONS**

See Known Interfering Substances and Stability and Storage Instructions under SAMPLE.

Keep plasma samples at room temperature to avoid low temperature activation.

It is recommended to respect the incubation time of samples with Reagent A.

Improper handling of samples may result in partial activation of coagulation factors which would cause erroneous results. Lupus anticoagulant may affect factor activity determination. A new calibration is required for each reagent batch and for each instrument used.

**PERFORMANCE**

a) **Reproducibility:** was determined using different samples

(in series and day-to-day). The following results were obtained:

**Intra-assay precision:**

Level	S.D.	C.V.
22.9 sec	0.460 sec	2.01%
21.2 sec	0.148 sec	0.70%
24.3 sec	0.579 sec	2.39%

**Inter-assay precision:**

Level	S.D.	C.V.
31.1 sec	0.802 sec	2.58%
32.9 sec	0.643 sec	1.95%
31.4 sec	1.513 sec	4.82%

**b) Measuring range:** 6-150%

**PARAMETERS FOR AUTOMATIC ANALYZERS**

Refer to the specific applications for each instrument.

**WIENER LAB PROVIDES**

5 x 1 ml (Code 1705015).

**REFERENCES**

- Nicolaes G et al (2002) Factor V and Thrombotic Disease. Description of a Janus-Faced Protein, ATVB, 22: 530-538.
- Duckers et al (2009) Advances in understanding the bleeding diathesis in factor V deficiency, BJH, 146: 17-26.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5<sup>th</sup> ed. CLSI: H21-A5.
- Triplett DA, (1981) New Methods in Coagulation, Crit Rev Clin Lab Sci. 1981;15 (1):25-84.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4<sup>th</sup> ed. 2001.

**SYMBOLS**

The following symbols are used in the packaging for Wiener lab. diagnostic reagents kits.

 This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

 Authorized representative in the European Community

 "In vitro" diagnostic medical device

 Contains sufficient for <n> tests

 Use by

 Temperature limitation (store at)

 Do not freeze

 Biological risks

 Volume after reconstitution

 Contents

 Batch code

 Manufactured by:

 Harmful

 Corrosive / Caustic

 Irritant

 Consult instructions for use

 Calibrator

 Control

 Positive Control

 Negative Control

 Catalog number

 Wiener Laboratorios S.A.I.C.  
 Riobamba 2944  
 2000 - Rosario - Argentina  
<http://www.wiener-lab.com.ar>  
 Tech. Dir.: Viviana E. Cétola  
 Biochemist  
 A.N.M.A.T. Registered Product  
 PM-1102-112

 **Wiener lab.**

2000 Rosario - Argentina