Fructosamina

Control - 2 niveles

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PROVIDED REAGENTS

Control Level 1: solution containing fructosamine concentrations between 200-400 $\mu\text{mol/l}$ (1.7 and 3.5 mmol/l DMF, deoximorpholinofructose), with adequate preservatives. Values change from lot to lot.

Control Level 2: solution containing fructosamine concentrations between 400-800 $\mu\text{mol/l}$ (3.5 and 7 mmol/l DMF, deoximorpholinofructose), with adequate preservatives. Values change from lot to lot.

WARNINGS

Reagent are for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: are stable in refrigerator (2-10°C) until the expiration date shown on the box.

INSTABILITY OR DETERIORATION OF REAGENTS

Any variation in the control' organoleptic characteristics may indicate its deterioration.

REQUIRED MATERIAL (non-provided)

Refer to Wiener lab's **Fructosamina AA líquida** package insert.

PROCEDURE

Use in the same way as an unknown sample, according to the instructions included in Wiener lab's **Fructosamina AA líquida** kit. The designated values for each level are indicated in the table below.

PROCEDURE LIMITATIONS

Refer to PROCEDURE LIMITATIONS in Wiener lab's **Fructosamina AA líquida** package insert.

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Wiener Laboratorios S.A.I.C.