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# HbA1c

## Calibrator

For the calibration of the HbA1c determination using immunoturbidimetric methods

### USES

**HbA1c Calibrator Turbitest AA** is designed for the calibration of the HbA1c determination using the **HbA1c v.2 Turbitest AA** kit from Wiener lab. Refer to the table of assigned values, since they are lot-specific.

### PROVIDED REAGENTS

**Calibrator:** lyophilized glycated and non-glycated human hemoglobin.

### NON-PROVIDED REAGENTS

Bidistilled or deionized water

### INSTRUCTIONS FOR USE

Open the vial by slowly removing the rubber cap to avoid loss of lyophilized material.

Add the deionized water quantity indicated on the vial label. Cap, let stand for 30 minutes, and then mix gently by inversion avoiding foam formation

The Calibrator does not require pretreatment with Wiener lab.'s **Reactivo Hemolizante**.

### WARNING

Reagents are for "in vitro" diagnostic use.

Do not ingest. Avoid contact with the skin. In case of exposure, proceed according to the instructions of the corresponding health authorities.

Use gloves and adequate eye protection.

The calibrator has been prepared from non-reactive material for hepatitis B (HBsAg) surface antigen, antibodies against hepatitis C virus (HCV) and antibodies against human immunodeficiency virus (HIV). However, since no method can rule out the potential risk of infection with absolute certainty, the controls and all the samples should be handled as if they were infective material.

Use the reagents according to the working procedures for clinical laboratories.

All reagents and samples should be discarded according to current regulations.

### STABILITY AND STORAGE INSTRUCTIONS

**Calibrator:** stable in refrigerator (2-10°C) until the expiration date stated on the box. After reconstitution it is stable for 8 hours at room temperature (15-25°C), for 2 days refrigerated (2-10°C) or for 3 months frozen (at -20°C). Aliquot and freeze the calibrator immediately after reconstitution. It must only be thawed once.

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Wiener Laboratorios S.A.I.C.