HbA1c

enzymatic - Control

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HbA1c

enzymatic - Control

For the precision control of HbA1c determination

USE

HbA1c enzymatic Control is designed to be used for the quality control of **HbA1c enzymatic** kit.

Refer to the table of assigned values for the constituents, since they are lot-specific.

PROVIDED REAGENTS

Control level 1: lyophilized human hemoglobin with HbA1c level within the range of normal values or within the clinically expected values for controlled diabetic patients.

Control level 2: lyophilized human hemoglobin with HbA1c level within the range of pathological values.

NON-PROVIDED REAGENTS

Demineralized water.

INSTRUCTIONS FOR USE

Open the vial by slowly removing the rubber cap to avoid loss of lyophilized material.

Add the demineralized water quantity indicated on the vial label.

Cap and mix by gentle inversion several times for a period of 30 minutes, avoiding foam formation.

Controls require pretreatment using **HbA1c enzymatic Lysis Buffer**.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Do not ingest. Avoid contact with skin and eyes.

Controls were prepared from non-reactive material for HBsAg, HCV and HIV. However, controls should be handled as if they were infective samples.

Use the reagents according to the working procedures for clinical laboratories.

All reagents and samples should be discarded according to current regulations.

STABILITY AND STORAGE INSTRUCTIONS

Lyophilized controls are stable at 2-10°C until the expiration date stated on the box. After reconstitution they are stable for up to 14 days at 2-10°C. After use, keep bottles tightly capped and stored at 2-10°C.

INSTABILITY OR DETERIORATION OF REAGENTS

Reagents must be limpid. Discard if turbidity is observed.

PROCEDURE

Controls must be used in the same way as an unknown sample, according to the instructions for use included in the **HbA1c enzymatic** kit. Remove from refrigerator, reconstitute and mix carefully because vigorously agitation can create bubbles. Transfer aliquots of the required volume to the cuvette of the instrument to be used. Cap and store at 2-10°C.

%HbA1c values may be expressed according to IFCC or DCCT/NGSP standardization using the following calculation methods:

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Wiener Laboratorios S.A.I.C.