HDL Colesterol

FT

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HDL Cholesterol

FT

Precipitating reagent (phosphotungstic acid) for HDL-cholesterol's separation in serum or plasma

SUMMARY

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins.

Phospholipids, free cholesterol and proteins constitute the outer surface of the lipoprotein particle, while the inner core contains mostly esterified cholesterol and triglycerides. Lipoproteins solubilize and transport cholesterol into the bloodstream.

Lipoproteins classification is based on its density which depends on the proportion of protein and lipid involved. Classes are: chylomicron, very-low-density lipoproteins (VLDL), low-density lipoproteins (LDL) and high-density lipoproteins (HDL). Numerous clinical studies have shown that the different lipoprotein classes have diverse effects on coronary heart disease risk.

The main role of HDL in lipid metabolism is the uptake and transport of cholesterol from peripheral tissues to the liver through a process known as reverse cholesterol transport (cardioprotective mechanism). Low HDL-cholesterol is associated with a high risk of coronary heart disease. Therefore, the determination of HDL-cholesterol is a useful tool for identification high-risk patients.

PRINCIPLE

High-density lipoproteins (HDL) are selectively separated by precipitating the low and very-low-density lipoproteins (LDL and VLDL) with phosphotungstic acid in magnesium ions presence. HDL remain in the supernatant. Cholesterol associated to HDL is determined by the enzymatic system cholesterol-oxidase/peroxidase with Trinder colorimetry (Phenol/4-Amino-phenazone).

PROVIDED REAGENTS

A. Reagent A (Precipitating Reagent): 0.44 mmol/l phosphotungstic acid and 20 mmol/l magnesium chloride solution.

NON-PROVIDED REAGENTS

Wiener lab.'s **Colestat enzimático**, **Colestat enzimático AA** or **Colestat enzimático AA líquida**.

INSTRUCTIONS FOR USE

Reagent A: ready to use.

WARNINGS

The Reagent is for "in vitro" diagnostic use. Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

The Provided Reagent is stable at room temperature (< 25°C) until the expiration date stated on the box.

INSTABILITY OR DETERIORATION OF REAGENTS

Any trace of bacterial contamination may indicate reagent deterioration.

SAMPLE

Serum or plasma

- a) Collection:** obtain the sample in the usual way.
- b) Additives:** if plasma is used, collect only with heparin.
- c) Known interfering substances:** anticoagulants other than heparin and bilirubin over 50 mg/l interfere. See Young, D.S. in References for effect of drugs on the present method.
- d) Stability and storage instructions:** separate serum within one hour from extraction. The Lipid Research Clinics recommend to refrigerate the sample until the test is performed. The conservation of the samples at room temperature alters their lipoprotein composition, even before 24 hours. Some authors mention that sample stability is maintained for up to 3 days at 4°C and it is extended when they are frozen. However, there is a great variability between different samples, therefore, it is recommended to keep them at 2-10°C and performe the test within 24 hours.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer or photocolormeter.
- Micropipettes and pipettes for measuring the stated volumes
- Kahn tubes.
- Spectrophotometric cuvettes.
- Water bath at 37°C.
- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 505 nm in spectrophotometer or 490-530 nm in photocolormeter with green filter.
- Reaction temperature: 37°C
- Reaction time: 35 minutes
- Sample volume: 200 ul
- Reagent A volume: 500 ul
- Supernatant volume: 200 ul
- Working Reagent volume of **Colestat enzimático**, **Colestat enzimático AA** or **Colestat enzimático AA líquida**: 2 ml
- Final reaction volume: 2.2 ml

PROCEDURE

In a Kahn tube adds 200 ul sample and 500 ul Reagent A. Homogenize by shaking (not by inversion) for 20 seconds and let repose for 10 minutes at room temperature. Centrifuge 15 minutes at 3000 rpm or 2 minutes at 12,000 rpm. Use clean supernatant as sample.

In three tubes labeled B (Blank), S (Standard) and U (Unknown) place:

	B	S	U
Supernatant	-	-	200 ul
Standard	-	20 ul	-
Working Reagent	2 ml	2 ml	2 ml

When Working Reagent of **Colesterol enzimático AA** or **Colestat enzimático AA líquida** is used, mix and incubate for up to 5 minutes at 37°C or 15 minutes at 37°C when Working Reagent of **Colestat enzimático** is used. Then, take the tubes out of the bath and bring them to room temperature.

Read absorbance at 505 nm in spectrophotometer or at 490-530 nm in photocolimeter with green filter, setting the instrument to zero with the Blank.

STABILITY OF FINAL REACTION

The reaction color is stable for up to 2 hours, thus the absorbance must be measure within that period.

CALCULATIONS

$$\text{HDL Cholesterol (g/l)} = U \times f = \frac{0.762}{S}$$

$$0.762 = 2 \text{ (g/l)} \times \frac{FV_E}{V_S} \times \frac{RV_E}{RV_{Std}} \times \frac{V_{Std}}{V_E}$$

FV_E = extract final volume = 0.7 ml

V_S = processed sample volume = 0.2 ml

RV_E = reaction volume with extract = 2.2 ml

RV_{Std} = reaction volume with Standard = 2.02 ml

V_{Std} = Standard volume in the reaction = 0.020 ml

V_E = extract volume in the reaction = 0.2 ml

If Reagent volumes other than 2 ml are used, the 0.762 factor varies and must be calculated again, replacing VR_E and VR_{Std} in the formula.

REFERENCE VALUES

The National Cholesterol Education Program (NCEP) Expert Panel provided the following reference values:

0.40 - 0.60 g/l

It is recommended that each laboratory establishes its own reference values. However, values over 0.40 g/l are recommended and the ones over 0.60 g/l have been regarded as protective. On the contrary, HDL-cholesterol values below 0.40 g/l is considered as risk indicator of coronary heart disease.

PROCEDURE LIMITATIONS

- See Known interfering substances under SAMPLE.
- The accuracy and precision of the determination mainly depend on following the precipitation conditions.
- Samples with triglycerides above 10 g/l may difficult the precipitation yielding turbid supernatants or a lipoprotein layer floating on the surface. In that case, dilute the sample 1:2 with saline solution and repeat the precipitation. The obtained result must be multiplied by 2.

PERFORMANCE

a) Reproducibility: processing replicates of one sample on the same day, the following results were obtained:

Level	S.D.	C.V.
0.32 g/l	± 0.011 g/l	3.4 %
0.68 g/l	± 0.024 g/l	3.5 %

b) Linearity: the reaction is linear up to 5 g/l.

c) Detection limit: depends on the spectrophotometer used. For a 0.001 O.D. reading, the minimum visible change of concentration will be of approximately 0.0063 g/l.

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Wiener Laboratorios S.A.I.C.