# Homocysteine

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(Cód. 1009931)\*

## BIBLIOGRAFIA

- McCully KS. - Am. J. Pathol. 56:111-128; 1969.
- McCully KS. - Ann. Clin. Lab. Sci. 23:447-493; 1993.
- Cramer DA. - Laboratory Medicine - 29:7; 1998.
- Rosenquist TH, Ratashak SA, Selhub J. - Proc. Natl. Acad. Sci. USA 93:15227-15232; 1996.
- Young, DS. Effects of drugs on clinical laboratory tests, 5<sup>th</sup> ed., AACCPress, Washington, DC, p 3-441; 2000.
- Ueland PM, Refsum H, Stabler SP, et al, - Clin. Chem. 39:1764-79; 1993.
- Probst R, Brandl R, Blümke M, Neumeier D. - Clin. Chem. 44:1567-9; 1998.
- CLSI EP9 A2: Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Ed., CLSI, 2002.
- CLSI EP15 A2: User verification of Performance for Precision and Trueness; Approved Guideline - Second Ed., CLSI, 2005.
- CLSI EP6 A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, CLSI, 2003.
- Silvia Persichilli, Bruno Zappacosta, et al. - Clin. Chem. Lab. Med. 46/12: 1786-1788; 2008.
- John M Scott et al. - Clin. Chem. 50/1:3-32; 2004.

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(Cód. 1009931)\*

#### REFERÊNCIA

- McCully KS. - Am. J. Pathol. 56:111-128; 1969.
- McCully KS. - Ann. Clin. Lab. Sci. 23:447-493; 1993.
- Cramer DA. - Laboratory Medicine - 29:7; 1998.
- Rosenquist TH, Ratashak SA, Selhub J. - Proc. Natl. Acad. Sci. USA 93:15227-15232; 1996.
- Young, DS. Effects of drugs on clinical laboratory tests, 5<sup>th</sup> ed., AACCPress, Washington, DC, p 3-441; 2000.
- Ueland PM, Refsum H, Stabler SP, et al. - Clin. Chem. 39:1764-79; 1993.
- Probst R, Brandl R, Blümke M, Neumeier D. - Clin. Chem. 44:1567-9; 1998.
- CLSI EP9 A2: Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Ed., CLSI, 2002.
- CLSI EP15 A2: User verification of Performance for Precision and Trueness; Approved Guideline - Second Ed., CLSI, 2005.
- CLSI EP6 A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, CLSI, 2003.
- Silvia Persichilli, Bruno Zappacosta, et al. - Clin. Chem. Lab. Med. 46/12: 1786-1788; 2008.
- John M Scott et al. - Clin. Chem. 50/1:3-32; 2004.

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# Homocysteine

For total homocysteine determination in serum and plasma


## SUMMARY

Homocysteine is a by-product of protein metabolism, specifically the intracellular demethylation of the essential amino acid, methionine. Homocysteine is exported into plasma where it circulates, mostly in its oxidized form, bound to plasma proteins.

Total homocysteine represents the sum of all homocysteine species found in serum or plasma (free plus protein bound). In small amounts, homocysteine is not harmful to the body or blood vessels, but when excess amounts accumulate in the blood stream, arterial vessels may be damaged and the resulting inflammation may eventually cause blockage of blood to the heart. Recent studies are yielding evidence that elevated blood levels of homocysteine have a predictive value for risk of coronary artery disease similar to that of elevated cholesterol levels. Recent evidence has also related elevated blood levels of homocysteine to the risk of miscarriages and birth defects.

## PRINCIPLE

Oxidized forms of homocysteine are reduced to free homocysteine which then reacts with serine catalyzed by cystathionine  $\beta$ -synthase (CBS) to form L-cystathionine. Cystathionine is cycled back to homocysteine by cystathionine  $\beta$ -lyase (CBL) forming in the process pyruvate and ammonia. This cycling sequence increases the sensitivity of the method by about 10-fold. Pyruvate is detected using lactate dehydrogenase (LDH) with NADH. The rate of conversion of NADH to NAD is directly proportional to the amount of homocysteine in the sample or calibrator.

## PROVIDED REAGENTS

**A. Reagent A:** > 35 kU/l lactate dehydrogenase (rabbit muscle), serine 0.76 mmol/l, NADH 0.47 mmol/l, Tris buffer pH 8.5.

**B. Reagent B:** > 20 kU/l cystathionine  $\beta$ -synthase (microbial), > 10 kU/l cystathionine  $\beta$ -lyase (microbial), pH 7.5.

## NON-PROVIDED REAGENTS

- Saline Solution (0.9%)
- Wiener lab's **Homocysteine Calibrator**

## INSTRUCTIONS FOR USE

**Provided Reagents:** ready to use.

## WARNINGS

The reagents are for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable at 2-10°C until the expiration date stated on the box.

## INSTABILITY OR DETERIORATION OF REAGENTS

Reagents are limpid solutions. Discard in case of turbidity.

## SAMPLE

Serum or plasma

**a) Collection:** obtain serum in the usual way.

**b) Additives:** in case plasma is used as sample, use EDTA or heparin for its collection.

**c) Known interfering substances:** no interference is observed by bilirubin up to 20 mg/dl, hemoglobin up to 500 mg/dl and triglycerides up to 500 mg/dl.

Refer to Young bibliography for effects of drugs on this method.

**d) Stability and storage instructions:** homocysteine in serum or plasma is stable for up to 4 days at 20-25°C, for up to 4 weeks at 2-10°C and for up to 4 years at -20°C.

## REQUIRED MATERIAL (non-provided)

- Volumetric material for measuring stated volumes.
- Automated analyzer.

## PROCEDURE

(Automated analyzer)

The general testing procedure for **Homocysteine** in an automated analyzer is detailed below. For programming instructions check the user's manual of the automated analyzer in use.

<b>Sample or Calibrator</b>	10 $\mu$ l
<b>Reagent A</b>	150 $\mu$ l

Incubate for 300 seconds at 37°C.

**Reagent B**

15 µl

Incubate for 60 seconds at 37°C. Measure initial absorbance at 340 nm ( $A_1$ ). At exactly 300 seconds measured (using stopwatch) the absorbance again ( $A_2$ ). To obtain homocysteine result in µmol/l, multiply the absorbance difference ( $\Delta A = A_2 - A_1$ ) by the factor.

**CALIBRATION**

Wiener lab's **Homocysteine Calibrator** is recommended. Calibrator should be processed as a sample and the corresponding factor is calculated based on it. Insert the calibrator's concentration value every time a lot is changed.

**QUALITY CONTROL METHOD**

For each determination, process 2 levels of a quality control material (**Homocysteine Control**) with known homocysteine concentrations.

Each laboratory should establish its own acceptable control range.

**REFERENCE VALUES**

Pregnant: < 10 mmol/l

Children < 15 years: < 10 mmol/l

Adults 15-65 years: < 15 mmol/l

Adults > 65 years: < 20 mmol/l

Consider these limits as guidelines only.

Homocysteine levels can vary with age, gender, geographical area and genetic factors, therefore it is important for laboratories to establish their own reference ranges.

**PROCEDURE LIMITATIONS**

See Known Interfering Substances under SAMPLE.

Do not pipette by mouth.

It is recommended to use Wiener lab's **Homocysteine Control** as quality control material. The use of controls from other manufacturers may yield different values for certain ranges because they are dependent from the method or system used.

**PERFORMANCE**

**a) Reproducibility:** the following coefficients of variations were obtained as an estimation of the intra-assay (CV<sub>i</sub>) and total (CV<sub>t</sub>) precision studies, according to the guidelines contained in CLSI EP15-A document:

**Intraassay precision (n = 20)**

Level	S.D.	C.V.
7.06 mmol/l	± 0.40 mmol/l	5.6%
12.38 mmol/l	± 0.38 mmol/l	3.1%
29.16 mmol/l	± 0.84 mmol/l	2.9%

**Total precision (n = 20)**

Level	S.D.	C.V.
7.06 mmol/l	± 0.38 mmol/l	5.4%
12.38 mmol/l	± 0.41 mmol/l	3.3%
29.16 mmol/l	± 1.39 mmol/l	4.8%

**b) Linearity:** reacción es linear up to homocysteine concentration of 60 mmol/l.

**c) Detection Limit:** 0.274 mmol/l.

**PARAMETERS FOR AUTOMATIC ANALYZERS**

For programming instructions check the user's manual of the automated analyzer in use.

**WIENER LAB. PROVIDES**

55 ml: - 1 x 50 ml Reagent A  
- 1 x 5 ml Reagent B  
(Cat. N° 1009365)

55 ml: - 1 x 50 ml Reagent A  
- 1 x 5 ml Reagent B  
(Cat. N° 1009623)

55 ml: - 1 x 50 ml Reagent A  
- 1 x 5 ml Reagent B  
(Cat. N° 1009931)\*

**REFERENCES**

- McCully KS. - Am. J. Pathol. 56:111-128; 1969.
- McCully KS. - Ann. Clin. Lab. Sci. 23:447-493; 1993.
- Cramer DA. - Laboratory Medicine - 29:7; 1998.
- Rosenquist TH, Ratashak SA, Selhub J. - Proc. Natl. Acad. Sci. USA 93:15227-15232; 1996.
- Young, DS. Effects of drugs on clinical laboratory tests, 5<sup>th</sup> ed., AACC Press, Washington, DC, p 3-441; 2000.
- Ueland PM, Refsum H, Stabler SP, et al, - Clin. Chem. 39:1764-79; 1993.
- Probst R, Brandl R, Blümke M, Neumeier D. - Clin. Chem. 44:1567-9; 1998.
- CLSI EP9 A2: Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Ed., CLSI, 2002.
- CLSI EP15 A2: User verification of Performance for Precision and Trueness; Approved Guideline - Second Ed., CLSI, 2005.
- CLSI EP6 A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, CLSI, 2003.
- Silvia Persichilli, Bruno Zappacosta, et al. - Clin. Chem. Lab. Med. 46/12: 1786-1788; 2008.
- John M Scott et al. - Clin. Chem. 50/1:3-32; 2004.

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(Nr kat. 1009931)

#### ŹRÓDŁA

- McCully KS. - Am. J. Pathol. 56:111-128; 1969.
- McCully KS. - Ann. Clin. Lab. Sci. 23:447-493; 1993.
- Cramer DA. - Laboratory Medicine - 29:7; 1998.
- Rosenquist TH, Ratashak SA, Selhub J. - Proc. Natl. Acad. Sci. USA 93:15227-15232; 1996.
- Young, DS. Effects of drugs on clinical laboratory tests, 5th ed., AACCC Press, Washington, DC, p 3-441; 2000.
- Ueland PM, Refsum H, Stabler SP, et al, - Clin. Chem. 39:1764-79; 1993.
- Probst R, Brandl R, Blümke M, Neumeier D. - Clin. Chem. 44:1567-9; 1998.
- CLSI EP9 A2: Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Ed., CLSI, 2002.
- CLSI EP15 A2: User verification of Performance for Precision and Trueness; Approved Guideline - Second Ed., CLSI, 2005.
- CLSI EP6 A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, CLSI, 2003.
- Silvia Persichilli, Bruno Zappacosta, et al. - Clin. Chem. Lab. Med. 46/12: 1786-1788; 2008.
- John M Scott et al. - Clin. Chem. 50/1:3-32; 2004.

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Wiener Laboratorios S.A.I.C.