# ISE TAR *plus* iones

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For ion determination using autoanalyzer's I.S.E. module

## PRINCIPLE

The automated determination of electrolyte concentration, such as Na<sup>+</sup>, K<sup>+</sup> and Cl<sup>-</sup> in serum or plasma, is based on a potentiometric method using ion-selective electrodes. Such electrodes consist of membranes that are permeable to a single ionic species. According to Nernst's equation, the potential difference established in the membrane interface and the sample in the solution, is directly proportional to the logarithm of the activity or ionic concentration:

$$E = E^{\circ} + RT/nF \times \log A$$

Where E is the potential difference and A is the ionic activity in the solution.

## PROVIDED REAGENTS

**Buffer:** 0.1 mol/l borate buffer, pH 8.3.

**Cleaning Solution:** 0.1 mol/l sodium hydroxide.

**Reference Solution:** 0.1 mol/l borate buffer with Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup> and CO<sub>3</sub>H<sup>-</sup> in low concentrations, pH 8.3.

**High Calibrator:** Na<sup>+</sup>, K<sup>+</sup> and Cl<sup>-</sup> solution in high concentrations. See the concentration stated on the label.

**Low Calibrator:** Na<sup>+</sup>, K<sup>+</sup> and Cl<sup>-</sup> solution in low concentrations. See the concentration stated on the label.

### Deproteinizing solution:

**Catalyzer:** envelope containing pepsin for a final concentration of 8 g/l.

**Solution:** 0.075 N hydrochloric acid.

## INSTRUCTIONS FOR USE

**Buffer, Reference Solution, Cleaning Solution, High and Low Calibrators:** ready to use.

**Deproteinizing solution:** add the content of a Catalyzer envelope to a Solution bottle. Mix well.

## WARNINGS

- Reagents are for "in vitro" diagnostic use.
- Do not expose Calibrators to air for more than 15 minutes. Calibrators may be deteriorated when they are exposed to air. Therefore, it is recommended to pour them into the cuvette where the samples are placed and calibrate immediately.
- Samples may be deteriorated in contact with air.
- Buffer should be immediately capped after use since it may be completely deteriorated by air contact in 3 days.
- The Cleaning Solution and the Solution (Deproteinizing solution) are corrosive. H315+H320: Causes skin and eye irritation. H314 Causes severe skin burns and eye damage. P262 Do not get in eyes, on skin, or on clothing. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and

easy to do. Continue rinsing. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P280 Wear protective gloves/protective clothing/eye protection/face protection.

- Use reagents according to the working procedures for clinical laboratories.
- Reagents and samples should be discarded according to the local regulations in force.

## STABILITY AND STORAGE INSTRUCTIONS

**Provided Reagents:** stable at room temperature (< 25°C) until their expiration date stated on the box.

**Deproteinizing solution:** once prepared it is stable for up to one month at 2-10°C

## REQUIRED MATERIAL (non-provided)

Autoanalyzer with ISE module.

## PROCEDURE

For use, refer to the autoanalyzer's user manual to be used.

## INSTRUMENT CALIBRATION

For calibration use the provided High and Low Calibrators. It should be performed every 4 hours or at least twice a day. It should also be performed in the following situations:

- when the quality control is outside the acceptable range;
  - every time the Buffer or the Reference solution are replaced.
- Calibration variations should be within the limits of the calibration curve showed in the following table:

Analyte	Calibration curve limit (mv)
Sodium	50 to 70
Potassium	40 to 70
Chloride	-60 to -35

## REFERENCE VALUES

Serum	Reference Range	ISE Linearity
Na <sup>+</sup>	134-149 mEq/l	20-400 mEq/l
K <sup>+</sup>	3.6-5.5 mEq/l	2.0-200 mEq/l
Cl <sup>-</sup>	94-112 mEq/l	40-400 mEq/l

## QUALITY CONTROL METHOD

Good laboratory practices require every day testing for at

least two control's levels. If control values are within the established higher and lower limits, run is considered under control and results may be reported. If control is not found within acceptable ranges, repeat the run using a new control and fresh reagents. Recalibrate if the results are still outside the established limits.

### WIENER LAB. PROVIDES

- Cat. Nr. 1009279:
- 5 x 10 ml Buffer
- 8 x 10 ml Reference Solution
- 5 x 10 ml Cleaning Solution
- 1 x 125 ml Solution
- 1 Catalyzer envelope
- 1 x 50 ml High Calibrator
- 1 x 50 ml Low Calibrator

### SYMBOLS EXPLANATION

<b>Sol</b>   <b>Referencia</b>	<b>Calibr</b>   <b>Alto</b>
Reference Solution	High Calibrator
<b>Sol</b>   <b>Limpieza</b>	<b>Calibr</b>   <b>Bajo</b>
Cleaning Solution	Low Calibrator
<b>Solución</b>	<b>Buffer</b>
Solution	Buffer
<b>Catalizador</b>	
Catalyzer	

### SYMBOLS

The following symbols are used in the packaging for Wiener lab. diagnostic reagents kits.

 This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

-  Authorized representative in the European Community
-  "In vitro" diagnostic medical device
-  Contains sufficient for <n> tests
-  Use by
-  Temperature limitation (store at)
-  Do not freeze
-  Biological risks
-  Volume after reconstitution
-  Contents
-  Batch code
-  Manufactured by:
-  Harmful
-  Corrosive / Caustic
-  Irritant
-  Consult instructions for use
-  Calibrator
-  Control
-  Positive Control
-  Negative Control
-  Catalog number

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Wiener Laboratorios S.A.I.C.