λ -Light Chain

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(Cód. 1009965)*

BIBLIOGRAFIA

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.

- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.

- Lievens M. - Medical and technical usefulness of measurement of kappa and lambda immunoglobulin light chains in serum with an M-component. - J. Clin. Chem. Clin. Biochem 27:519-523, 1989.

- Hafner G et al. - Effects of standardization with the new international reference preparation for proteins in human serum on method comparability and reference values. - Clin. Lab. 41:743-748, 1995.

- Jones RG et al. - Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: I. Detection. - Clin. Chem. 37:1917-1921, 1991.

- Jones RG et al. - Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: II. Classification by use of Computer-based algorithms. - Clin. Chem. 37:1922-1926, 1991.

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(Cód.1009965)*

REFERÊNCIAS

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACCPress, 5th ed., 2000.
- Lievens M. - Medical and technical usefulness of measurement of kappa and lambda immunoglobulin light chains in serum with an M-component. - J. Clin. Chem. Clin. Biochem 27:519-523, 1989.
- Hafner G et al. - Effects of standardization with the new international reference preparation for proteins in human serum on method comparability and reference values. - Clin. Lab. 41:743-748, 1995.
- Jones RG et al. - Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: I. Detection. - Clin. Chem. 37:1917-1921, 1991.
- Jones RG et al. - Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: II. Classification by use of Computer-based algorithms. - Clin. Chem. 37:1922-1926, 1991.

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λ -Light Chain

Immunoturbidimetric method for quantitative determination of λ -light chain in serum

SUMMARY

The polyclonal immunoglobulins exhibit kappa and lambda light chains in a constant ratio of 2:1, while monoclonal immunoglobulins show only a single light chain type. Increased production of monoclonal immunoglobulins or free monoclonal light chains produces a ratio kappa / lambda outside the reference range, indicating the existence of a monoclonal gammopathy.

The determination of light chains is useful for diagnosis and monitoring of multiple myeloma, lymphocytic neoplasms, Waldenstrom's macroglobulinemia and connective tissue diseases such as rheumatoid arthritis or lupus erythematosus (SLE).

PRINCIPLE

The λ -light chain reacts to the specific antibody forming insoluble immune complexes. The turbidity caused by these immune complexes is proportional to λ -light chain concentration in the sample and may be spectrophotometrically measured.

PROVIDED REAGENTS

A. Reagent A: phosphate buffer, pH 7.4.

B. Reagent B: polyclonal antibodies anti-human λ -light chain (goat) in phosphate buffer, pH 7.4.

NON-PROVIDED REAGENTS

- Saline solution
- Wiener lab.'s **Calibrador Proteínas nivel alto Turbitest AA**

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

The reagents are for "in vitro" diagnostic use. All patient samples should be handled as though capable of transmitting infectious diseases. Use the reagents according to the working procedures for clinical laboratories. The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date stated on the box. Do not freeze.

SAMPLE

Serum

a) Collection: obtain in the usual way.

b) Additives: not required.

c) Known interfering substances: do not use hemolyzed, lipemic or contaminated samples. Before testing, particles in samples should be removed by centrifugation.

No interferences have been observed with hemoglobin up to 1100 mg/dl, triglycerides up to 1200 mg/dl, bilirubin up to 20 mg/dl and rheumatoid factor up to 300 U/l.

See Young, D.S. in References for effect of drugs on the present method.

d) Stability and storage instructions: sample should be preferably fresh. In case it cannot be processed immediately, sample can be kept for up to 48 hours at 2-10°C or for longer period store at -20°C.

REQUIRED MATERIAL (non-provided)

- Spectrophotometer.
- Square spectrophotometric cuvettes.
- Micropipettes and pipettes for measuring the stated volumes.
- Kahn or hemolysis tubes.
- Stopwatch.

ASSAY CONDITIONS

- Wavelength: 340 nm
- Reaction temperature: room temperature (25°C). Temperature control is not critical, it can range between 22 and 30°C.
- Reaction time: 15 minutes
- Sample volume: 20 μ l
- Final reaction volume: 1120 μ l

Sample and reagents volumes may be proportionally changed without affecting the calculation factors.

PROCEDURE

CALIBRATION CURVE

In Kahn tubes dilute the **Calibrador Proteínas nivel alto** with saline solution 1:10, 1:20, 1:40, 1:80 and 1:160, using saline solution as the zero point.

Diluted Calibrador Proteínas	20 μ l
-------------------------------------	------------

Reagent A	1000 μ l
------------------	--------------

Homogenize and measure the absorbance of each dilution at 340 nm (OD₃₄₀), setting the instrument to zero with distilled water. Then, add:

Reagent B	100 μ l
------------------	-------------

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD_2), setting the instrument to zero with distilled water.

Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each Calibrador Proteínas dilution, including the zero point.

Draw on graph paper the ΔA absorbance differences based on the Calibrador Proteínas concentration in mg/dl (g/l).

SAMPLES PROCEDURE

Dilute the Samples 1:10 with saline solution.

Diluted Sample	20 ul
Reagent A	1000 ul

Homogenize and measure the absorbance at 340 nm (OD_1), setting the instrument to zero with distilled water. Then add:

Reagent B	100 ul
------------------	--------

Mix and incubate 15 minutes at room temperature. Measure the absorbance at 340 nm (OD_2), setting the instrument to zero with distilled water.

Intra-assay precision

Level	S.D.	C.V.
55.5 mg/dl	± 1.2 mg/dl	2.1 %
140.0 mg/dl	± 2.0 mg/dl	1.4 %
333.2 mg/dl	± 5.3 mg/dl	1.6 %

b) Detection limit: 25 mg/dl.

c) Measuring range: 25 - 400 mg/dl.

d) Prozone effect: not noted until 8600 mg/dl λ -light chain.

WIENER LAB. PROVIDES

65 ml: - 1 x 60 ml Reagent A
- 1 x 5 ml Reagent B
(Cat. N° 1009358)

65 ml: - 1 x 60 ml Reagent A
- 1 x 5 ml Reagent B
(Cat. N° 1009655)

65 ml: - 1 x 60 ml Reagent A
- 1 x 5 ml Reagent B
(Cat. N° 1009965)*

REFERENCES

- Dati, F et al. - Proteins-Laboratory testing and clinical use, 2005.
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 5th ed., 2000.
- Lievens M. - Medical and technical usefulness of measurement of kappa and lambda immunoglobulin light chains in serum with an M-component. - J. Clin. Chem. Clin. Biochem 27:519-523, 1989.
- Hafner G et al. - Effects of standardization with the new international reference preparation for proteins in human serum on method comparability and reference values. - Clin. Lab. 41:743-748, 1995.
- Jones RG et al. - Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: I. Detection. - Clin. Chem. 37:1917-1921, 1991.
- Jones RG et al. - Use of Immunoglobulin Heavy-chain and Light-chain measurement in a multicenter trial to investigate Monoclonal components: II. Classification by use of Computer-based algorithms. - Clin. Chem. 37:1922-1926, 1991.

CALCULATIONS

Calculate the absorbance difference ($\Delta A = OD_2 - OD_1$) for each sample tested. Interpolate this ΔA in the calibration curve to determine the concentration in mg/dl (g/l) corresponding to the sample under study. Samples with an absorbance above that of the **Calibrador Proteínas nivel alto** must be diluted with saline solution and processed again. Multiply the obtained result by the dilution factor.

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Control should be processed in the same manner as samples.

REFERENCE VALUES

93 - 242 mg/dl (0.93 - 2.42 g/l)

Ratio κ/λ : 1.17 - 2.93

Each laboratory should set its own reference values.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

It is recommended to perform a complete recalibration when changing reagent lot or when suggested by Quality Control.

Avoid contamination to preserve the integrity of the reagents. Only use thoroughly clean and dry micropipettes for measurement.

PERFORMANCE

a) Reproducibility: replicates of samples containing different λ -light chain levels were assayed and the following results were obtained:

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Wiener Laboratorios S.A.I.C.