Lithium

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1 x 20 mL Reagente A (Cód. 1009972)

REFERÊNCIAS

- Chapoteau, E., Czech, B.P., Zazulak, W. and Kumar, A. First Colorimetric Assay of Lithium in Serum. (1992), Clin. Chem. 9: 1654-1657.
- Moyer, T.P. and Pippenger, C.E. Therapeutic Drug Monitoring. Tietz Textbook of Clinical Chemistry. Burtis, C. A. and Ashwood, E.R. eds, (Second ed.) W.B. Saunders Company.
- Amdisen, A. Serum Lithium Determinations for Clinical Use (1967), Scand. J. Clin. Lab Invest. 20: 104-108.
- Wachtel, M.S. et al, Creation and Verification of Reference Intervals (1995), Laboratory Medicine 26: 593-597.
- Linder, M.W. and Keck, P.E Jr. Standards of Laboratory Practice: Antidepressant drug monitoring. (1998) Clin. Chem. 44:5 1073-1084.
- "Use of Anticoagulants in Diagnostic Investigations". (Jan 2002) WHO Publication WHO/DIL/LAB 99.1 Rev 2.

Lithium

For the quantitative determination of lithium in serum and plasma

SUMMARY

Lithium, administered as lithium carbonate, is used for the treatment of the manic phase of affective disturbances, mania and manic-depressive illness. It acts by increasing the recapture of neurotransmitters, thus reducing their concentration in the neuronal junction. This produces a sedative effect on the central nervous system. It is completely absorbed by the gastrointestinal tract and maximum serum levels occur 2 to 4 hours after an oral dose. The serum half-life is 48 to 72 hours and is eliminated through the kidneys.

Reduction of renal function may prolong the clearance time. The first symptoms of intoxication include apathy, slowness, drowsiness, lethargy, speech difficulties, irregular tremors, myoclonic contractions, muscle weakness and ataxia. Levels greater than 1.5 mmol/L (12 hours after dose) indicate a significant risk of poisoning. Lithium serum concentrations are essentially performed to ensure compliance and avoid toxicity.

PRINCIPLE

The method is based on a complex formation among lithium ions with a lithium-specific chromoionophore in alkaline solution. The concentration of lithium in the sample is proportional to the absorbance variation by the complex formation and is measured at 510/480 nm.

PROVIDED REAGENTS

A. Reagent A: sodium hydroxide (0.5 mol/L), EDTA (50 µmol/L), chromoionophore (15 µmol/L), preservative and detergent.

NON-PROVIDED REAGENTS

Wiener lab.'s **Calibrador A plus**.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Do not pipette with your mouth.

The reagent is corrosive. H315+H320: Causes skin and eye irritation. H314 Causes severe skin burns and eye damage. P262 Do not get in eyes, on skin, or on clothing. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P280 Wear protective gloves/protective clothing/eye protection/face protection.

Use the reagent keeping the usual working precautions in

the clinical analysis laboratory.

All reagents and samples must be discarded according to current local regulations.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable at 2-10°C until the expiration date stated on the box, protected from light.

SAMPLE

Serum or plasma

a) Collection: collect the sample in the usual way. In order to standardize the test and evaluate the therapy properly, it is recommended to collect the sample 12 hours after receiving the last dose of lithium.

b) Additives: in case of using plasma, the use of sodium heparin or dipotassium EDTA is recommended as an anticoagulant.

c) Known interfering substances:

Criteria: recovery within $\pm 10\%$ of the value of the initial values to therapeutic connections.

No interference is observed by conjugated bilirubin up to 30 mg/dL, unconjugated bilirubin up to 30 mg/dL, hemoglobin up to 1000 mg/dL, triglycerides up to 1500 mg/dL. No significant interference with other cations was found at physiological concentrations for NaCl (140 mM), NH₄Cl (20 µmol/L), KCl (4 mM), CaCl₂ (2.5mM), MgCl₂ (1.0 mM), FeCl₃ (1 mg/L) and CuCl₂ (1.2 mM), ZnCl₂ (1 mM).

Refer to Young's bibliography for the effects of drugs in the present method.

d) Stability and storage instructions: serum is stable for 5 days at room temperature (20-25°C) tightly closed so as not to concentrate the sample by evaporation, for 30 days at 2-10°C or for 6 months frozen (-20°C).

REQUIRED MATERIAL (non-provided)

- Volumetric material to measure the indicated volumes
- Automated analyzer

PROCEDURE

(automated analyzers)

A general procedure for **Lithium** in an automated analyzer is detailed below. When the technique is implemented for a particular analyzer, follow its work instructions.

In a cuvette maintained at the chosen temperature, place:

Sample or Calibrator

4 µL

Reagent A200 μ L

Incubate for 10 minutes at 37°C. Absorption reading at 510 nm as λ_1 and 450 or 480 nm as λ_2 against water blank.

CALCULATIONS

[Lithium] mmol/L = (Δ Abs) x Factor

Wiener lab. analyzers automatically calculate the lithium concentration of each sample.

CALIBRATION

Calibrador A plus is processed in the same way as the samples and the corresponding factor is calculated from it. It is recommended to use two-point calibration after changing reagent lot and when required by quality control.

QUALITY CONTROL METHOD

With each determination, process 2 levels of a quality control material (**Standatrol S-E 2 niveles** from Wiener lab.) with known concentrations of lithium.

REFERENCE VALUES

Therapeutic range: 0.60 - 1.20 mmol/L

Alarm range: 1.20 - 1.50 mmol/L

Risk of intoxication: > 1.50 mmol/L

Severe poisoning: > 2.50 mmol/L

It is recommended that each laboratory establish its own reference values.

Toxicity can occur even when serum levels are within the recommended therapeutic range. For this reason it is important to determine the optimal range for each patient individually.

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

PERFORMANCE

a) Accuracy: was evaluated according to CLSI protocol EP5.

In this study two samples were used with different levels of lithium tested with 2 daily runs in duplicate for 20 days.

Repeatability (n = 80)

Level	S.D.	C.V.
0.62 mmol/L	0.016 mmol/L	2.5%
1.63 mmol/L	0.015 mmol/L	0.9%

Intra-laboratory precision (n = 80)

Level	S.D.	C.V.
0.62 mmol/L	0.018 mmol/L	2.8%
1.63 mmol/L	0.038 mmol/L	2.4%

b) Analytical sensitivity: the minimum detectable concentration of lithium was determined at 0.05 mmol/L.

c) Measuring range: 0.05 - 3.00 mmol/L

d) Correlation: lithium value was determined in 100 serum samples with a range between 0.05 and 3.00 mmol/L with **Lithium** from Wiener lab. and a commercial kit based on the

same principle, obtaining the following correlation coefficient: $r = 0.9974$; slope $b = 0.9813 x$ and intersection $a = 0.0154$

PARAMETERS FOR AUTOMATED ANALYZERS

For programming instructions refer to the user manual of the analyzer in use.

WIENER LAB. PROVIDES

1 x 20 mL Reagent A (Cat. N° 1999810)

1 x 20 mL Reagent A (Cat. N° 1009705)

1 x 20 mL Reagent A (Cat. N° 1009972)

REFERENCES

- Chapoteau, E., Czech, B.P., Zazulak, W. and Kumar, A. First Colorimetric Assay of Lithium in Serum. (1992), Clin. Chem. 9: 1654-1657.
- Moyer, T.P. and Pippenger, C.E. Therapeutic Drug Monitoring. Tietz Textbook of Clinical Chemistry. Burtis, C. A. and Ashwood, E.R. eds, (Second ed.) W.B. Saunders Company.
- Amdisen, A. Serum Lithium Determinations for Clinical Use (1967), Scand. J. Clin. Lab Invest. 20: 104-108.
- Wachtel, M.S. et al, Creation and Verification of Reference Intervals (1995), Laboratory Medicine 26: 593-597.
- Linder, M.W. and Keck, P.E Jr. Standards of Laboratory Practice: Antidepressant drug monitoring. (1998) Clin. Chem. 44:5 1073-1084.
- "Use of Anticoagulants in Diagnostic Investigations". (Jan 2002) WHO Publication WHO/DIL/LAB 99.1 Rev 2.

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Wiener Laboratorios S.A.I.C.