PCR

Control N

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For the quality control of C-reactive protein (CRP)
determination by immunoturbidimetric method

USES

PCR Control N Turbitest AA has been designed for the quality control of Wiener lab's **CRP hs Turbitest AA**. Check the value stated on the label, since it is lot-specific.

PROVIDED REAGENTS

PCR Control N: human C-reactive protein solution containing CRP levels within normal values or within the limit of normal/pathological values.

INSTRUCTIONS FOR USE

Provided Reagent: ready to use.

Before testing, gently invert the vial several times avoiding foam formation. Dispense the Control volume required to be used. Do not return the remaining material into the original vial.

This Control has been specifically designed to be used as PCR ultrasensible Turbitest AA kit's control.

WARNINGS

The reagent is for "in vitro" diagnostic use.

Avoid ingestion and contact with the skin.

The Control has been prepared from non-reactive material to HBsAg, HCV and HIV. However, Control and all samples should be handled as infectious material.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

The unopened Control is stable at 2-10°C until the expiration date stated on the box.

Once opened, it is stable for up to 2 weeks at 2-10°C.

Tightly cap the vial and store it at 2-10°C. Protect from light.

When stored at -20°C, the control is stable for up to 3 months.

Once opened it is recommended to store at -20°C in fractions.

Do not repeatedly freezing and thawing.

INSTABILITY OR DETERIORATION OF REAGENTS

Turbidity or microbial growth indicates contamination and deterioration of the reagent. Discard in such case.

PROCEDURE LIMITATIONS

Avoid reagent contamination with any action such as: transferring the content by inversion, placing any device into the vial, etc.

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Wiener Laboratorios S.A.I.C.