



Tensioactive

Tensioactive agent to be used in automated analyzers

USES

Tensioactive is employed in automated analyzers in order to obtain a better rinsing of the hydraulic circuit fluids, avoiding bubble formation that can have an impact on the precision of results.

PROVIDED REAGENTS

Tensioactive: Tween 20® 1% solution. Non-ionic tensioactive, biodegradable.

NON PROVIDED REAGENTS

Demineralized or distilled water.

INSTRUCTIONS FOR USE

Add **Tensioactive** to demineralized or distilled water according to the specifications given for each analyzer.

WARNINGS

The reagent is for "in vitro" diagnostic use.
Use the reagents according to the working procedures for clinical laboratories.
All reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Tensioactive: stable at 2-10°C until the expiration date stated on the box.

INSTABILITY AND DETERIORATION OF REAGENTS


Evident signs of turbidity or presence of particles in suspension indicate reagent deterioration. Discard in such case. Avoid external contamination.

WIENER LAB. PROVIDES

- 1 x 200 ml (Cat. 1999719).


SYMBOLS

The following symbols are used in packaging for Wiener lab. diagnostic reagents kits.

 This product fulfills the requirements of the European Directive 98/79 EC for "in vitro" diagnostic medical devices

 Authorized representative in the European Community

 "In vitro" diagnostic medical device


 Contains sufficient for <n> tests

 Use by


 Temperature limitation (store at)

 Do not freeze

 Biological risks


 Volume after reconstitution

 Contents


 Batch code

 Manufactured by:


 Harmful

 Corrosive / Caustic

 Irritant

 Consult instructions for use


 Calibrator

 Control

 Positive Control

 Negative Control

 Catalog number

 Wiener Laboratorios S.A.I.C.
Riobamba 2944
2000 - Rosario - Argentina
<http://www.wiener-lab.com.ar>
Dir. Téc.: Viviana E. Cétola
Bioquímica
Producto Autorizado A.N.M.A.T.
PM-1102-46



Wiener lab.

2000 Rosario - Argentina