WL 19 Con

AA

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WL 19 Con

AA

For quality control in hematology analyzers

USES

WL 19 Con AA is designed to be used in the quality control of specific hematology parameters in Wiener lab. Counter 19 and Wiener lab. Counter 19 CP analyzers. Please refer to the table of assayed values, since they are lot-specific.

PROVIDED REAGENTS

Control: suspension containing three different levels of human erythrocytes, artificial platelet components and artificial leukocytes with adequate preservatives.

INSTRUCTIONS FOR USE

Control: ready to use.

Before testing, gently invert the vial several times. Ensure that the content of the vial is completely suspended. Do not shake.

WARNINGS

Reagents are for "in vitro" diagnostic use.

The Controls have been prepared from non-reactive material for Hepatitis B surface antigen (HBsAg), antibodies to Hepatitis C virus (HCV) and antibodies to the human immunodeficiency virus (HIV). Nonetheless, the controls and samples should be handled as capable of transmitting infection.

Use the reagents according to the working procedures for clinical laboratories.

All reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

The control is stable at 2-10°C until the expiration date stated on the box. Do not freeze.

Once opened, it is stable for up to 14 days at 2-10°C. After use, close vials tightly and keep them at 2-10°C. Avoid contamination.

INSTABILITY AND DETERIORATION OF REAGENTS

In the non-homogenized Control vials the supernatant is reddish and slightly turbid. Once homogenized, the normal color is red to dark red. Other organoleptic feature variation (like the dark red supernatant) or unacceptable results may indicate deterioration. Discard in such case.

PROCEDURE

1- Verify that the analyzer is adequately prepared, as stated in the operation manual.

2- Remove the vial from refrigerator and bring to room temperature (15-30°C) for 15 minutes before mixing. Do not use a mechanical mixer.

3- Roll the vial with the palm of the hands in horizontal position for 20-30 seconds. Occasionally invert the vial. Mix well, but do not shake.

4- Repeat step 3 until a complete suspension is obtained. Vials stored for long periods of time may require extra mixing.

5- Gently invert the vial 8-10 times.

6- Make sure the vial content is completely suspended. Otherwise, repeat the steps 3 to 5 until a total suspension is obtained.

7- Clean residual material from the cap and tube rim with a lint-free tissue.

8- Cap immediately after use. Do not keep at room temperature for more than 30 minutes. Store at 2-10°C.

EXPECTED RESULTS

Verify that the lot number on the package matches the lot number on the instructions for use. Assay values are determined on well-maintained, properly calibrated instruments using the manufacturer's recommended reagents.

Each laboratory should establish its own values and acceptable ranges and periodically reevaluate the mean.

PROCEDURE LIMITATIONS

A leukocyte manual differential analysis cannot be performed with this control.

Incomplete mixing of a vial prior to use invalidates both the sample withdrawn and any remaining material in the vial.

Do not use after the expiration date stated on the label.

Reagent differences, instrument maintenance, operation technique and calibration, may contribute to intra-laboratory variation.

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Wiener Laboratorios S.A.I.C.