

CK-MB

Control - 3 niveles

LOT 1905314920

LOT Nivel 1 312650

LOT Nivel 2 312670

LOT Nivel 3 312680

VALORES ASIGNADOS

VALORES ASIGNADOS DE CK-MB (U/l)* - BT 3000 PLUS / WIENER LAB, CB SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	28.5	20.0	37.1	40.7	30.5	50.9	87.2	65.4	109.0
CK-MB LÍQ.	29.5	20.7	38.4	43.3	32.5	54.1	90.1	67.6	112.6

VALORES ASIGNADOS DE CK-MB (U/l)* - KONELAB SERIES / WIENER LAB. CT SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	32.1	22.5	41.7	49.5	37.1	61.9	100	75	125
CK-MB LÍQ.	31.9	22.3	41.5	47.4	35.6	59.3	97.8	73.4	122.3

VALORES ASIGNADOS DE CK-MB (U/l)* - METROLAB 2300 PLUS / WIENER LAB. CM SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	26.8	18.8	34.8	40.2	30.2	50.3	82.2	61.7	102.8
CK-MB LÍQ.	32.7	22.9	42.5	47.8	35.9	59.8	91.9	68.9	114.9

VALORES ASIGNADOS DE CK-MB (U/l)* - TARGA BT 3000									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	28.5	20.0	37.1	40.7	30.5	50.9	87.2	65.4	109.0
CK-MB LÍQ.	29.5	20.7	38.4	43.3	32.5	54.1	90.1	67.6	112.6

VALORES ASIGNADOS DE CK-MB (U/l)* - WIENER LAB. CMD 800 SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB LÍQ.	29.2	20.4	38.0	45.9	34.4	57.4	91.9	68.9	114.9

* a 37°C

** Valores obtenidos empleando técnica con calibración / Valores obtidos utilizando técnica com calibração / Values obtained using calibration technique / Wartości uzyskane stosując techniki z kalibracją

Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
kit	kit	kit	zestaw
nivel 2	nível 2	level 2	poziom 2
nivel 3	nível 3	level 3	poziom 3
valores asignados	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango acceptable	faixa aceitável	acceptable range	pożądany zakres

