

CMD Urine

Quality Control

LOT

1901292320

Rev 000

LOT

Nivel 1 291400

Fecha 20190114

LOT

Nivel 2 291410

VALORES ASIGNADOS

Na ⁺	Valor Medio meq/l	Rango meq/l	
Nivel 1	74.2	63.0	85.3
Nivel 2	143	122	165

K ⁺	Valor Medio meq/l	Rango meq/l	
Nivel 1	20.2	17.2	23.2
Nivel 2	45.3	38.5	52.1

Cl ⁻	Valor Medio meq/l	Rango meq/l	
Nivel 1	92.7	78.8	106.5
Nivel 2	183	155	210

Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
nivel 1	nível 1	level 1	poziom 1
nivel 2	nível 2	level 2	poziom 2
valores	valores	values	wartości
asignados	estabelecidos	assigned	zakres
valor medio	valor médio	mean	średnia wartość
rango	faixa	range	zakres
fecha	data	date	data

