

# Control Plasma

*Normal*

LOT 1904306570

Exp: 2021/04

## Reference Values

Reagent	Manual Technique		Automated Technique*	
Prothrombin Time (sec)	10.1	16.9	10.1	16.9
APTT (sec)	24.4	40.6	-	-
APTT ellagic (sec)	20.9	34.8	20.4	34.0
Fibrinogen (mg/dl)	260	390	260	390
Thrombin Time (sec)				
4.0 UNIH/ml	11.0	16.4	11.0	16.4
2.7 UNIH/ml	14.2	21.4	14.2	21.4

(\*) Coa DATA coagulation instruments

### Dutch Diagnostics BV

Hekkehorst 3 - 7207 BS Zutphen

The Netherlands

Phone: +31 575 512 300

Fax: +31 575 510 164

E-mail: [info@dutch-diagnostics.com](mailto:info@dutch-diagnostics.com)