

# Control Plasma

## *Pathologic*

LOT 1904306580

Exp: 2021/05

### Reference Values

Reagent	Manual Technique		Automated Technique*	
Prothrombin Time (sec)	26.5	44.1	26.5	44.1
APTT (sec)	47.3	78.8	-	-
APTT ellagic (sec)	41.6	69.4	41.3	68.9
Fibrinogen (mg/dl)	116	174	116	174
Thrombin Time (sec)				
4.0 UNIH/ml	12.9	19.3	12.9	19.3
2.7 UNIH/ml	16.2	24.2	16.2	24.2

(\*) Coa DATA coagulation instruments

#### **Dutch Diagnostics BV**

Hekkehorst 3 - 7207 BS Zutphen

The Netherlands

Phone: +31 575 512 300

Fax: +31 575 510 164

E-mail: [info@dutch-diagnostics.com](mailto:info@dutch-diagnostics.com)