

Fructosamina

LOT

1903303500

LOT

Nivel 1

303090

LOT

Nivel 2

303080

VALORES ASIGNADOS - TECNICA MANUAL / ANALIZADORES						
TECNICA/ ANALIZADOR	NIVEL 1			NIVEL 2		
	VALOR MEDIO	RANGO		VALOR MEDIO	RANGO	
Konelab Series / Wiener lab, CT Series	296 umol/L	237	355	467 umol/L	373	560
	2.6 mmol/L (DMF)	2.1	3.1	4.1 mmol/L (DMF)	3.2	4.9
Targa - BT 3000	282 umol/L	226	338	490 umol/L	392	588
	2.5 mmol/L (DMF)	2.0	2.9	4.3 mmol/L (DMF)	3.4	5.1
Metrolab 2300 plus / Wiener lab. CM Series	275 umol/L	220	330	467 umol/L	374	561
	2.4 mmol/L (DMF)	1.9	2.9	4.1 mmol/L (DMF)	3.3	4.9
BT 3000 plus / Wiener lab. CB Series	282 umol/L	226	338	490 umol/L	392	588
	2.5 mmol/L (DMF)	2.0	2.9	4.3 mmol/L (DMF)	3.4	5.1
CMD 800 Series	296 umol/L	237	355	476 umol/L	381	571
	2.6 mmol/L (DMF)	2.1	3.1	4.1 mmol/L (DMF)	3.3	5.0
Técnica manual	299 umol/L	239	359	568 umol/L	454	681
	2.6 mmol/L (DMF)	2.1	3.1	4.9 mmol/L (DMF)	4.0	5.9

Glosario/Glossário/Glossary/Glosariusz

Español

Português

English

Polski

nivel
valores asignados
rango
analizadores

nível
valores estabelecidos
faixa
analisadores

level
assigned values
range
analyzers

poziom
zakres wartości
zakres
analizator

**Wiener lab.**