

Fructosamina

LOT 2004370360

Control - 2 niveles

LOT Nivel 1 367830

LOT Nivel 2 367840

VALORES ASIGNADOS - TECNICA MANUAL / ANALIZADORES						
TECNICA/ ANALIZADOR	NIVEL 1			NIVEL 2		
	VALOR MEDIO	RANGO		VALOR MEDIO	RANGO	
Konelab Series / Wiener lab, CT Series	213 umol/L	171	256	403 umol/L	322	483
	1.9 mmol/L (DMF)	1.5	2.2	3.5 mmol/L (DMF)	2.8	4.2
Targa - BT 3000	198 umol/L	158	237	410 umol/L	328	491
	1.7 mmol/L (DMF)	1.4	2.1	3.6 mmol/L (DMF)	2.8	4.3
Metrolab 2300 plus / Wiener lab. CM Series	196 umol/L	157	235	404 umol/L	323	485
	1.7 mmol/L (DMF)	1.4	2.0	3.5 mmol/L (DMF)	2.8	4.2
BT 3000 plus / Wiener lab. CB Series	198 umol/L	158	237	410 umol/L	328	491
	1.7 mmol/L (DMF)	1.4	2.1	3.6 mmol/L (DMF)	2.8	4.3
CMD 800 Series	206 umol/L	165	247	393 umol/L	314	472
	1.8 mmol/L (DMF)	1.4	2.2	3.4 mmol/L (DMF)	2.7	4.1
Técnica manual	172 umol/L	137	206	440 umol/L	352	528
	1.5 mmol/L (DMF)	1.2	1.8	3.8 mmol/L (DMF)	3.1	4.6

Glosario/Glossário/Glossary/Glosariusz

Español

nivel
valores asignados
rango
analizadores

Português

nível
valores estabelecidos
faixa
analísadores

English

level
assigned values
range
analyzers

Polski

poziom
zakres wartości
zakres
analizator



Wiener lab.