

Proti U/LCR

Control - 2 niveles

LOT 2006379530

LOT Nivel 1 377990

LOT Nivel 2 378000

VALORES ASIGNADOS PARA TECNICA MANUAL Y ANALIZADORES						
TECNICA/ ANALIZADOR	NIVEL 1			NIVEL 2		
	VALOR MEDIO (mg/dL)	RANGO (mg/dL)		VALOR MEDIO (mg/dL)	RANGO (mg/dL)	
Konelab Series / Wiener lab, CT Series	28.8	21.6	35.9	96.6	77.3	115.9
Metrolab 2300 plus / Wiener lab. CM Series	30.9	23.2	38.7	93.5	74.8	112.1
Targa BT 3000	31.9	24.0	39.9	94.3	70.7	117.9
BT 3000 plus / Wiener lab. CB Series	31.9	24.0	39.9	94.3	70.7	117.9
Wiener lab. CMD Series	29.0	21.8	36.3	94.6	75.7	113.5
Técnica Manual	36.3	27.2	45.4	97.1	72.8	121.4

Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
nivel	nível	level	poziom
fecha	data	date	data
valores asignados	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango	faixa	range	zakres
técnica manual	técnica manual	manual technique	metody manualne
analizadores	analísadores	analyzers	analizator

