

Proti U/LCR

Control - 2 niveles

LOT 2007384170

LOT Nivel 1 384050

LOT Nivel 2 384060

VALORES ASIGNADOS PARA TECNICA MANUAL Y ANALIZADORES						
TECNICA/ ANALIZADOR	NIVEL 1			NIVEL 2		
	VALOR MEDIO (mg/dL)	RANGO (mg/dL)		VALOR MEDIO (mg/dL)	RANGO (mg/dL)	
Konelab Series / Wiener lab, CT Series	32.1	24.1	40.2	99.3	79.5	119.2
Metrolab 2300 plus / Wiener lab. CM Series	34.6	26.0	43.3	100	80	120
Targa BT 3000	35.0	26.3	43.8	98.4	73.8	123.0
BT 3000 plus / Wiener lab. CB Series	35.0	26.3	43.8	98.4	73.8	123.0
Wiener lab. CMD Series	32.7	24.5	40.9	101	81	121
Técnica Manual	35.6	26.7	44.5	102	77	128

Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
nivel	nível	level	poziom
fecha	data	date	data
valores asignados	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango	faixa	range	zakres
técnica manual	técnica manual	manual technique	metody manualne
analizadores	analisadores	analyzers	analizator

