

CK-MB

Control - 3 niveles

LOT 1911349190

LOT Nivel 1 345870

LOT Nivel 2 345880

LOT Nivel 3 345890

VALORES ASIGNADOS

VALORES ASIGNADOS DE CK-MB (U/l)* - BT 3000 PLUS / WIENER LAB, CB SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	28.5	20.0	37.1	40.6	30.5	50.8	81.3	61.0	101.6
CK-MB LÍQ.	31.7	22.2	41.2	43.7	32.8	54.6	87.5	65.6	109.4

VALORES ASIGNADOS DE CK-MB (U/l)* - KONELAB SERIES / WIENER LAB. CT SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	33.0	23.1	42.9	46.7	35.0	58.4	93.8	70.4	117.3
CK-MB LÍQ.	32.6	22.8	42.4	45.3	34.0	56.6	94.4	70.8	118.0

VALORES ASIGNADOS DE CK-MB (U/l)* - METROLAB 2300 PLUS / WIENER LAB. CM SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	27.4	19.2	35.6	39.5	29.6	49.4	77.5	58.1	96.9
CK-MB LÍQ.	39.6	27.7	51.5	53.8	40.4	67.3	97.9	73.4	122.4

VALORES ASIGNADOS DE CK-MB (U/l)* - TARGA BT 3000									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	28.5	20.0	37.1	40.6	30.5	50.8	81.3	61.0	101.6
CK-MB LÍQ.	31.7	22.2	41.2	43.7	32.8	54.6	87.5	65.6	109.4

VALORES ASIGNADOS DE CK-MB (U/l)* - WIENER LAB. CMD 800 SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB LÍQ.	28.6	20.0	37.2	40.9	30.7	51.1	84.3	63.2	105.4

* a 37°C

** Valores obtenidos empleando técnica con calibración / Valores obtidos utilizando técnica com calibração / Values obtained using calibration technique / Wartości uzyskane stosując techniki z kalibracją

Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
kit	kit	kit	zestaw
nivel 2	nível 2	level 2	poziom 2
nivel 3	nível 3	level 3	poziom 3
valores asignados	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango acceptable	faixa aceitável	acceptable range	pożądany zakres

