

# CK-MB

## Control - 3 niveles

LOT 2008389200

LOT Nivel 1 380750

LOT Nivel 2 380760

LOT Nivel 3 380770

### VALORES ASIGNADOS

VALORES ASIGNADOS DE CK-MB (U/l)* - BT 3000 PLUS / WIENER LAB, CB SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	29.1	20.4	37.8	47.1	35.3	58.9	85.4	64.1	106.8
CK-MB LÍQ.	30.3	21.2	39.4	48.5	36.4	60.6	89.4	67.0	111.7

VALORES ASIGNADOS DE CK-MB (U/l)* - KONELAB SERIES / WIENER LAB. CT SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	33.6	23.5	43.7	53.1	39.8	66.4	98.7	74.0	123.4
CK-MB LÍQ.	33.4	23.4	43.4	53.2	39.9	66.5	96.3	72.2	120.4

VALORES ASIGNADOS DE CK-MB (U/l)* - METROLAB 2300 PLUS / WIENER LAB. CM SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	29.5	20.7	38.4	46.3	34.7	57.9	84.5	63.4	105.6
CK-MB LÍQ.	36.1	25.3	46.9	54.0	40.5	67.5	92.6	69.5	115.8

VALORES ASIGNADOS DE CK-MB (U/l)* - TARGA BT 3000									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB DS	29.1	20.4	37.8	47.1	35.3	58.9	85.4	64.1	106.8
CK-MB LÍQ.	30.3	21.2	39.4	48.5	36.4	60.6	89.4	67.0	111.7

VALORES ASIGNADOS DE CK-MB (U/l)* - WIENER LAB. CMD 800 SERIES									
KIT	NIVEL 1			NIVEL 2			NIVEL 3		
	VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE		VALOR MEDIO	RANGO ACCEPTABLE	
CK-MB LÍQ.	29.3	20.5	38.1	46.6	35.0	58.3	88.2	66.2	110.3

\* a 37°C

\*\* Valores obtenidos empleando técnica con calibración / Valores obtidos utilizando técnica com calibração / Values obtained using calibration technique / Wartości uzyskane stosując techniki z kalibracją

### Glosario/Glossário/Glossary/Glosariusz

Español	Português	English	Polski
kit	kit	kit	zestaw
nivel 2	nível 2	level 2	poziom 2
nivel 3	nível 3	level 3	poziom 3
valores asignados	valores estabelecidos	assigned values	zakres wartości
valor medio	valor médio	mean	średnia wartość
rango acceptable	faixa aceitável	acceptable range	pożądany zakres

