

Control Plasma

Normal

LOT

2101410870

Reference Values

Reagent	Manual Technique		Automated Technique*	
Prothrombin Time (sec)	9.2	14.4	9.4	14.6
APTT (sec)	25.8	40.4	-	-
APTT ellagic (sec)	22.5	35.3	20.2	31.6
Fibrinogen (mg/dl)	256	384	259	389
Thrombin Time (sec)				
4.0 UNIH/ml	12.9	19.3	12.0	18.0
2.7 UNIH/ml	16.8	25.2	16.7	25.1

(*) Coa DATA coagulation instruments

Dutch Diagnostics BV

Hekkehorst 3 - 7207 BS Zutphen

The Netherlands

Phone: +31 575 512 300

Fax: +31 575 510 164

E-mail: info@dutch-diagnostics.com